# EIAgen PSA Equi Kit

**REF** LI4027C

 96 tests

**Adaltis Italia S.p.A.**  
Via Cristoni, 12  
40033 Casalecchio di Reno – (BO) Italy  
Tel. +39-051-6136511 – Fax + 39-051-575280  
www.adaltis.com

**FOR IN VITRO DIAGNOSTIC USE ONLY** 

Store at 2...8 °C

## SYMBOLS USED ON LABELS

<b>MTPLATE</b>	microplate
<b>CONJ HRP 20x</b>	HRP Conjugate Concentrate
<b>BIOTIN Anti-PSA</b>	Biotin Anti-PSA
<b>CAL0</b>	Calibrator 0
<b>CAL1</b>	Calibrator 1
<b>CAL2</b>	Calibrator 2
<b>CAL3</b>	Calibrator 3
<b>CAL4</b>	Calibrator 4
<b>CAL5</b>	Calibrator 5
<b>CONTROL1</b> <b>CONTROL2</b>	Controls
<b>SOLN TMB</b>	TMB substrate
<b>SOLN STOP</b>	Stop solution
<b>WASH BUF 25X</b>	washing buffer concentrate
<b>LOT</b>	Lot number
<b>REF</b>	Catalogue Code
	Expiry date (Use by...)
 8°C	Temperature limitation (store at 2...8°C)
	Number of tests
	Keep away from sunlight
	Manufactured by
	Attention, See Instructions For Use
<b>IVD</b>	In vitro diagnostic medical device (In vitro diagnostic use)
	Biological risks (see SECTION 7- Warnings and Precautions).

### 1.0 INTENDED USE

The EIAgen PSA Equi KIT is intended for the quantitative determination of Total Prostate Specific Antigen (PSA) in human serum.

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

PSA is a 32 kDa single chain glycoprotein serine protease with a chymotrypsin like specificity produced by the secretory epithelium of the prostate gland (1). PSA is normally secreted into the seminal fluid and plays a functional role in the cleavage of the seminal vesicle proteins and the liquefaction of the seminal coagulum (2).

Only low levels of PSA are normally present in the blood stream, and increasing serum concentrations indicate prostatic pathology, including benign prostatic hyperplasia (BPH) and cancer of the prostate (PCa). Determination of PSA is now widely used for detection and management of patients with prostatic cancer and considered as the superior serological marker for cancer of the prostate (3,8).

PSA has been shown to form stable complexes with different antiproteases and the dominating portion of PSA in patient serum occurs in complex with alfa1-antichymotrypsin (PSA-ACT) (4). However there are large variations in the relation between total PSA and PSA-ACT complex between different individuals.

Studies have also indicated that the proportion of free PSA is higher in benign prostatic disease as compared to prostatic cancer (5). The antibodies in the EIAgen PSA Equi kit have been carefully evaluated and selected to give the same molar response for free PSA as for the PSA-ACT complex (7). The EIAgen PSA Equi KIT thus gives a true "total" PSA value independent of the individual variations of free and ACT complexed PSA.

### 3.0 PRINCIPLE OF THE TEST

The EIAgen PSA Equi kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with biotinylated anti-PSA monoclonal antibody and horseradish peroxidase (HRP) labelled tracer in Streptavidin coated microtiter strips. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra - methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour (Optical Density OD) is proportional to the amount of PSA present in the samples. The colour intensity is determined in a microtiter plate photometer at 620 nm (or optionally at 405 nm after addition of Stop Solution). Calibration curves are constructed for each assay by plotting OD value versus the concentration for each calibrator. The PSA concentrations of patient samples are then read from the calibration curve.

### 4.0 CONTENT OF THE KIT

#### 4.1 MICROPLATE (code 30PIA)

1 microplate of 12 strips of 8 breakable wells. The wells are coated with streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

#### 4.2 CONJUGATE (code 30HRP)

1 vial containing 0,75 mL of Stock solution (20x concentrate) of HRP conjugated to Anti-PSA monoclonal antibody from mouse, approximately 20 µg/mL. Contains preservatives.

To be mixed with Biotin Anti- PSA prior to use.

**Preparation of Conjugate (HRP- Anti-PSA) working solution:** Prepare the required quantity of conjugate working solution by mixing 50 µL of conjugate 20x concentrate with 1 mL of Biotin Anti- PSA per strip (see table below).

N° of strips	CONJUGATE µl	Biotin Anti- PSA mL
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	300	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of the Conjugate working solution.

**Alternative:** Pour the content of the Conjugate 20x concentrate, into the vial of Biotin Anti- PSA and mix gently. Make sure that all of the Conjugate is transferred to the vial of Biotin Anti- PSA.

**STORAGE:** the Conjugate working solution is stable for 3 weeks at 2–8°C in a sealed container. Do not prepare more Conjugate working solution than will be used within this period and make sure that it is stored properly.

#### 4.3 BIOTIN Anti PSA (code 30BT) BIOTIN|Anti-PSA

1 vial containing 15 mL of Biotin labelled Anti-PSA monoclonal antibody from mouse, approximately 1 µg/mL. Contains phosphate buffered saline (pH 7.2), bovine serum albumin, bovine immunoglobulin, blocking agents, Tween 20, an inert blue dye and 0.01% methyl-isothiazolone (MIT) as preservative.  
Ready to use.

#### 4.4 CALIBRATORS

6 vials containing 0.75 mL of calibrator. The calibrators contain Human PSA in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye and 0.01% methyl-isothiazolone (MIT) as preservative.  
Ready for use.

Calibrator	symbol	code	concentration
Calibrator 0	<span style="border: 1px solid black; padding: 0 2px;">CAL0</span>	30CAL0	0 µg/L
Calibrator 1	<span style="border: 1px solid black; padding: 0 2px;">CAL1</span>	30CAL1	1 µg/L
Calibrator 2	<span style="border: 1px solid black; padding: 0 2px;">CAL2</span>	30CAL2	2 µg/L
Calibrator 3	<span style="border: 1px solid black; padding: 0 2px;">CAL3</span>	30CAL3	10 µg/L
Calibrator 4	<span style="border: 1px solid black; padding: 0 2px;">CAL4</span>	30CAL4	30 µg/L
Calibrator 5	<span style="border: 1px solid black; padding: 0 2px;">CAL5</span>	30CAL5	60 µg/L

#### Reference material

The 1st International Standard 96/670 may be used as a reference standard. Values for PSA Calibrators and Controls were assigned against a set of in-house reference standards whose values are traceable to the 1st International Standard.

#### 4.5 CONTROLS

Control	symbol	code	
Control 1	<span style="border: 1px solid black; padding: 0 2px;">CONTROL1</span>	30CTRL1	1 vial of 0.75 ml
Control 2	<span style="border: 1px solid black; padding: 0 2px;">CONTROL2</span>	30CTRL2	1 vial of 0.75 ml

Each control contains Human PSA in a Tris-HCl buffered salt solution containing bovine serum albumin, and 0.01% methylisothiazolone (MIT) as preservative.  
Ready for use.

#### 4.6 TMB-SUBSTRATE (code TMBB) SOLN|TMB

1 vial containing 12 mL of ready-to-use Hydrogen peroxide and 3,3',5,5' tetra-methylbenzidine (TMB) .  
Ready to use.

#### 4.7 STOP SOLUTION (code STOPB) SOLN|STOP

1 vial containing 15 mL of 0.12 M hydrochloric acid.  
Ready to use.

#### 4.8 WASHING BUFFER 25X (code TLAVB) WASH|BUF|25X

1 bottle containing 50 mL of wash buffer concentrate. Wash buffer concentrate contains a Tris-HCl buffered salt solution with Tween 20. contains German II as preservative.

**Preparation of Wash Buffer working solution:** Pour the 50 mL Washing Buffer Concentrate into a clean container and dilute 25-fold by adding 1200 mL of distilled or deionised water to give a Wash Buffer working solution.

**STORAGE:** the stability of diluted wash buffer is 2 weeks at 2-25°C in a sealed container.

#### 5.0 STORAGE AND STABILITY AFTER THE FIRST OPENING

- Store kit components at 2...8° C and do not use after the reported expiry date.
- After use, the microplate should be resealed, the bottle caps replaced and tightened and the kit stored at 2...8°C.
- The opened kit should be used within expiry date indicated on labels.

#### 6.0 MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- **Microtiter plate shaker:** Shaking should be medium to vigorous. Longitudinal shaking approximately 200 strokes/min, oscillations 700—900/min.
- **Microtiter plate wash device:** Automatic plate wash capable of performing 1, and 6 washing cycles, or semi manual microtiter plate washing device connected to vacuum pump or water-jet vacuum and a liquid trap for retaining aspirated liquid. The Nunc Immuno-8 manual

strip washer is recommended if an automatic microtiter plate wash is not used.

- **Microtiter plate photometer:** with a wavelength of 620 nm and/or 405 nm; and 450nm for ultra sensitive determination (see par 9.2 items B – J) and an OD range of 0 to 3.0.
- **Precision pipettes:** with disposable plastic tips for dispensing microlitre volumes. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential. Pipettes for dispensing millilitre volumes.
- **Distilled or deionised water:** For preparation of Wash Solution
- **Others:** reagent reservoirs for multichannel pipettes, paper towels, and timer

## 7.0 WARNINGS OR PRECAUTIONS

### 7.1 SAFETY PRECAUTIONS

- All reagents in this kit are for in vitro diagnostic use only
- Only experienced laboratory personnel should use this test and handling should be agreement with GLP.
- Operators should wear gloves and protective clothing when handling any patient sera or serum based products.
- Avoid contact with reagents containing hydrogen peroxide, hydrochloric acid and preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is total of all detergents.
- Materials used in the preparation of human source reagents have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.
- Follow all local, state and national regulations for disposal of all waste material. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

### 7.2 TECHNICAL PRECAUTIONS

#### A. Correct use of reagents

- The performance data represented here were obtained using specific reagents listed in this package insert. Do not use reagents from other manufacturers in the kits.
- Do not use reagents from other EIAgen kits with this kit. Do not mix reagents from different kit lots.
- Do not dilute or adulterate the kit reagents, unless directed by the kit protocol.
- The **TMB Substrate** is very sensitive for contamination. The TMB HRP Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded. For optimal stability of the TMB Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).

#### B. Correct pipetting procedure

- Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate solution.

#### C. Correct washing procedure

- Pre-washing of the strips before the assay is necessary
- Microplate washing is important. Improperly washed wells will give erroneous results  
WASHING PROCEDURE: A careful washing procedure of the strips is essential. Ensure that each well is filled up completely to the top edge and that the aspiration of the wells between and after the washing cycles is complete and the wells are dry. If there is liquid left in the wells, invert the plate and tap it carefully against absorbing paper.  
Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles prior to and after each incubation step. The aspiration/wash device should not be left standing with the Wash Solution for long periods as the needles may get clogged, giving poor liquid delivery and suction.

#### D. Adherence to assay procedure and specifications

- The test protocol must be followed strictly. Observe the indicated incubation times and temperature and the washing procedure. These are

critical steps. Do not allow the wells to dry between incubations. Please refer also to section 9.0\_ASSAY PROCEDURE.

- Include the positive and negative controls in every test run to monitor for reagent stability and correct assay performance. The obtained OD values for calibrators and concentration for controls have to be always compared to the ones reported in QC sheet. Please refer also to section 10.1\_VALIDITY OF THE ASSAY AND 10.2\_ QUALITY CONTROL.
- Do not use the kit to determine values outside the range indicated in the Instruction for Use. Please refer also to section 10.5\_MEASUREMENTS RANGE.

## 8.0 SPECIMEN COLLECTION AND STORAGE

The kit is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8°C for 24 hours. For longer periods stores samples at -20°C or below. Samples should not be stored in a self-defrosting freezer and not be thawed and refrozen before analysis. Allow frozen samples to thaw slowly, preferably at 2–8°C over night and then bring the samples to room temperature before analysis.

## 9.0 ASSAY PROCEDURE

### 9.1 REAGENT PREPARATION

- Bring all reagents to room temperature (20...25°C). The assay should only be performed at temperatures between 20–25°C to obtain accurate results.
- Prepare washing buffer working solution according to par. 4.8
- Prepare conjugate working solution according to par 4.2;
- Select sufficient microwells for the test. Transfer the required number of microtiter plate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant and reseal carefully). Before starting to pipette calibrators and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
- Perform each determination in duplicate for both calibrators and patient samples. A calibration curve should be run with each assay.

**NOTE:** The Calibrator 1 shall be used when ultra sensitive measurement is desired. Please refer to at par 9.2 items B and H Option for full instructions.

### 9.2 PIPETTING AND INCUBATION STEPS

- Wash** each strip once with the Wash Buffer working solution: Do not wash more strips than can be handled within 30 min.
- Pipette** 25 µL in duplicate of the Calibrators (CAL 0, 2, 3, 4, 5) , Controls (C 1, C 2) and patient specimens into the strip wells;  
**NOTE** For ultrasensitive determination of PSA in the low range (0-10 µg/L), CAL 1 (1 µg/L) should be included in the calibration curve. It is then necessary to exclude the two highest calibrators CAL 4 (30 µg/L) and CAL 5 ( 60 µg/L).
- Add** 100 µL of Conjugate working solution to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.
- Incubate** the plate for 1 hours (± 10 min) at room temperature (20-25°C) with constant shaking of the plate using a microtiter plate shaker.
- Aspirate** the incubation solution and **wash** each strip 6 times with the Wash Buffer working solution using the wash procedure described in section 7.2.C
- Add** 100 µL of TMB HRP-Substrate to each well using the same pipetting technique as in item C. The TMB Substrate should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.
- Incubate** for 30 min (± 5 min) at room temperature with constant shaking. Avoid exposure to direct sunlight.
- Immediately **read** the absorbance at 620 nm in a microtiter plate photometer.  
**Option 1:** If the laboratory does not have access to a microtiter plate photometer capable of reading at 620 nm, the OD can be determined as follows: Add 100 µL of Stop Solution. Mix and read the OD at 405 nm in a microtiter plate photometer within 15 minutes after addition of Stop Solution.  
**Option 2** For ultrasensitive determination of PSA in the low range 0-10 µg/L (see NOTE in item B), measurements of OD should be done according to Option 1 but reading at 450 nm.

## 10.0 CALCULATION OF RESULTS

### 10.1 VALIDITY OF THE ASSAY

Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated.

### 10.2 QUALITY CONTROL

Each laboratory may in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the precision of the assay.

If the results should fall out with the expected range, repeat the assay using freshly prepared controls. If the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall out with the acceptable ranges.

### 10.3 DATA REDUCTION: AUTOMATED METHOD

If a microtiter plate photometer reader with built-in data calculation program is used refer to the manual for the plate reader and create a program using the concentration stated on the labels of each Calibrator.

For automatic calculation of Total PSA results it is recommended to use either of the following methods:

- *Cubic spline* fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.
- *Spline smoothed* curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with *point-to-point* evaluation. Calibrator 0 should be included in the curve with the value 0 µg/L.
- *Quadratic curve* fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.

**NOTE:** 4-parametric or linear regression should not be used.

### 10.4 DATA REDUCTION: MANUAL METHOD

For manual evaluation, a calibration curve is constructed by plotting the OD values obtained for each PSA calibrator against the corresponding PSA concentration (in µg/L), see figure below. The unknown PSA concentrations can then be read from the calibration curve using the mean OD value of each patient specimen.

### 10.5 MEASUREMENT RANGE

The EIAgen PSA Equi kit measures concentrations between 0.1 and 60 µg/L. If PSA concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human male serum prior to analysis.

**NOTE:** The serum used for dilution should also be measured in order to determine the endogenous PSA concentration

If samples in an initial analysis give PSA levels higher than 60 µg/L the samples should be diluted 1/10 with normal male human serum and reanalyzed to obtain the accurate PSA concentration.

The PSA concentration of the undiluted sample is calculated as:

$$\text{Dilution } 1/10 = 10 \times ([\text{PSA}]_{\text{diluted sample}} - (0.9 \times [\text{PSA}]_{\text{Normal male serum}}))$$

### 11.0 EXPECTED VALUES

Healthy males are expected to have PSA values below 4 µg/L. However, age related variations have been shown to occur (6, 8). It is recommended that each laboratory establish its own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc.

### 12.0 LIMITATION OF THE PROCEDURE.

The level of PSA cannot be used as absolute evidence for the presence or absence of malignant disease. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients. The PSA test should not replace any established clinical examination.

The calibrators of the EIAgen PSA Equi kit should not be used for recovery studies of PSA. For recovery studies it is recommended to use a highly elevated patient sample.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

### 13.0 PERFORMANCE CHARACTERISTICS

#### 13.1 PRECISION

Total precision was determined according to NCCLS guideline EP5-A (9) using four levels of frozen pooled human serum containing added PSA and six

different EIAgen PSA Equi reagent combinations. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days.

Sample	Replicates	Mean µg/L	Within-run SD µg/L	Within run CV %	Between day SD µg/L	Between day CV %
PSA 1	80	1.42	0.04	2.7	0.03	2.2
PSA 2	80	5.92	0.13	2.2	0.06	1.0
PSA 3	80	14.2	0.35	2.5	0.12	0.8
PSA 4	80	39.2	0.60	1.5	0.60	1.5

### 13.2 DETECTION LIMIT

The detection limit of the EIAgen PSA equi kit is < 0.1 µg/L defined as the concentration corresponding to the mean of the OD values of the calibrator 0 plus 2 standard deviations according to formula:

$$(2 \times \text{SD CAL 0}) / (\text{OD CAL 2} - \text{OD CAL 0}) \times 2 \mu\text{g/L}$$

### 13.3 RECOVERY

Spiked serum samples were prepared by adding aliquots of samples with highly elevated PSA to normal male serum samples. The recovery of the antigen was within ± 10% of the expected values. NOTE: Recovery studies should not be performed using the kit calibrators.

### 13.4 HOOK EFFECT

No hook effect has been noticed with samples up to 23000 µg/L.

**NOTE:** In very high samples the colour of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low OD at 620 nm, and in extreme cases the OD may fall within the calibration curve range and noticed as a hook.

### 13.5 LINEARITY

Patient samples were diluted with normal male human serum and analyzed. The obtained values were within ± 10% of the expected values.

### 13.6 SPECIFICITY

The EIAgen PSA Equi kit is based on two mouse monoclonal antibodies, PSA10 and PSA66, directed against two distinct epitopes exposed both in the PSA-ACT complex and in free PSA. This antibody combination provides an equimolar response for free PSA and the PSA-ACT complex (7). The NCCLS guideline EP7-P (10) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant (± 10%) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

### 14.0 AUTOMATION

Application protocols for the proper automation (Personal Lab and Nexgen) on the Adaltis microtiterstrips analyzers are available upon request at Adaltis directly.

### 15.0 SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors.

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different (± 50%) from OD reported on QC	<ul style="list-style-type: none"> <li>- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>- incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation)</li> <li>- error in washing or in photometer reading (suggestion: check operating or settings of respective instruments)</li> <li>- contamination of Substrate (suggestion: use only disposable and clean plastic containers)</li> </ul>
Low reproducible results	<ul style="list-style-type: none"> <li>- not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>- error in washing or in reading (suggestion: check operating or settings of respective instruments)</li> <li>- contamination of Substrate (suggestion: use only disposable and clean plastic containers)</li> <li>- pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)</li> </ul>
no colourimetric reaction after addition of substrate	<ul style="list-style-type: none"> <li>- some reagent not pipetted</li> <li>- strong contamination of conjugate or Substrate</li> <li>- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)</li> </ul>
too low reaction (too low ODs)	<ul style="list-style-type: none"> <li>- incorrect conjugate dilution</li> <li>- incubation time too short, incubation temperature too low</li> </ul>
too high reaction (too high ODs)	<ul style="list-style-type: none"> <li>- incorrect conjugate dilution</li> <li>- incubation time too long, incubation temperature too high</li> <li>- water quality for wash buffer insufficient (low grade of deionization)</li> <li>- insufficient washing (conjugates not properly removed)</li> </ul>
unexplainable outliers	<ul style="list-style-type: none"> <li>- contamination of pipettes, tips or containers</li> <li>- inconstant and insufficient washing (conjugates not properly removed)</li> </ul>
too high within-run CV%	<ul style="list-style-type: none"> <li>- reagents and/or strips not pre-warmed to Room Temp. prior to use</li> <li>- plate washer is not washing correctly (suggestion: clean washer head)</li> </ul>
too high between-run CV%	<ul style="list-style-type: none"> <li>- incubation conditions not constant (time, temperature)</li> <li>- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)</li> <li>- person-related variation</li> </ul>

### 16.0 REFERENCES

1. Wang MC, Valenzuela LA, Murphy GP, Chu TM (1979). Purification of a human prostate specific antigen. Invest Urol 17: 159-163.
2. Lilja H. (1985). A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J Clin Invest 76: 1899-1903.
3. Oesterling JE (1991). Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urology 145: 907-923.
4. Lilja H., Christensson A., Dahlén U., Matikainen M-T, Nilsson O., Pettersson K., Lövgren T. (1991). Prostate-specific antigen in serum occurs predominantly in complex with α1-antichymotrypsin. Clin Chem 37: 1618-1625.
5. Christensson A., Björk T., Nilsson O., Dahlén U., Matikainen M-T., Cockett ATK, Abrahamsson PA, Lilja J. (1993). Serum prostate specific antigen complexed to α1-antichymotrypsin as an indicator of prostate cancer. J Urology 150: 100-105.
6. Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate. Specific antigen concentration: An important clinical observation. Urol Clin North Am 20: 671-80, 1993a.
7. Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms PSA. Br J Cancer 75(6): 789-797.
8. P Price C., Allard J., Davies G., Dawney A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem 2001; 38: 188-216.
9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)

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antigen complexed to a1-antichymotrypsin as an indicator of prostate cancer. *J Urology* 150: 100–105.

6. Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate. Specific antigen concentration: An important clinical observation. *Urol Clin North Am* 20: 671–80, 1993a.
7. Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms PSA. *Br J Cancer* 75(6): 789–797.
8. P Price C., Allard J., Davies G., Dawnay A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. Pre-and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. *Ann Clin Biochem* 2001; 38: 188–216.
9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)

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### 16.0 BIBLIOGRAFÍA

1. Wang MC, Valenzuela LA, Murphy GP, Chu TM (1979). Purification of a human prostate specific antigen. Invest Urol 17: 159-163.
2. Lilja H. (1985). A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J Clin Invest 76: 1899-1903.
3. Oesterling JE (1991). Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urology 145: 907-923.
4. Lilja H., Christensson A., Dahlén U., Matikainen M-T, Nilsson O., Pettersson K., Lövgren T. (1991). Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha 1$ -antichymotrypsin. Clin Chem 37: 1618-1625.
5. Christensson A., Björk T., Nilsson O., Dahlén U., Matikainen M-T., Cockett ATK, Abrahamsson PA, Lilja (1993). Serum prostate specific antigen complexed to  $\alpha 1$ -antichymotrypsin as an indicator of prostate cancer. J Urology 150: 100-105.
6. Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate. Specific antigen concentration: An important clinical observation. Urol Clin North Am 20: 671-80, 1993a.
7. Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms PSA. Br J Cancer 75(6): 789-797.
8. P Price C., Allard J., Davies G., Dawney A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem 2001; 38: 188-216.
9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)

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#### 16.0 RIFERIMENTI

- Wang MC, Valenzuela LA, Murphy GP, Chu TM (1979). Purification of a human prostate specific antigen. Invest Urol 17: 159-163.
- Lilja H. (1985). A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J Clin Invest 76: 1899-1903.
- Oesterling JE (1991). Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urology 145: 907-923.
- Lilja H., Christensson A., Dahlén U., Matikainen M-T, Nilsson O., Pettersson K., Lövgren T. (1991). Prostate-specific antigen in serum occurs predominantly in complex with α1-antichymotrypsin. Clin Chem 37: 1618-1625.
- Christensson A., Björk T., Nilsson O., Dahlén U., Matikainen M-T., Cockett ATK, Abrahamsson PA, Lilja H. (1993). Serum prostate specific antigen complexed to α1-antichymotrypsin as an indicator of prostate cancer. J Urology 150: 100-105.
- Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate. Specific antigen concentration: An important clinical observation. Urol Clin North Am 20: 671-80, 1993a.
- Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms of PSA. Br J Cancer 75(6): 789-797.
- P Price C., Allard J., Davies G., Dawnay A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. Pre-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem 2001; 38: 188-216.
- National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. *Ann Clin Biochem* 2001; 38: 188–216.

9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986)