# ImmuGlo™ Anti-Adrenal Cortex Antibody (AcAb) Test

IVD

PRODUCT INSERT

Code 1171      48 Determinations

## INTENDED USE

An indirect immunofluorescence antibody test for the detection and semi-quantitation of anti-adrenal cortex antibodies in human serum.

## SUMMARY AND EXPLANATION

Addison's disease is an endocrine disorder that occurs in all age groups. The disease is characterized by weight loss, muscle weakness, fatigue, low blood pressure and skin pigmentation. Most cases of Addison's disease are caused by the destruction of the adrenal cortex. About 70% of cases of Addison's disease are autoimmune mediated, while tuberculosis accounts for 20% and 10% are due to less common causes such as chronic infections, amyloidosis, metastatic neoplasia, and surgical removal of the adrenal gland<sup>1</sup>. Major clinical manifestations of Addison's disease (anorexia, abdominal pain, wasting, apathy, weakness, fasting hypoglycemia, diminished ability to conserve sodium and excrete free water, hyponatremia, increased production of ACTH, and  $\beta$ 2 lipoprotein) are attributable to deficiencies of cortisol and aldosterone<sup>2</sup>.

The diagnosis of Addison's disease is first made by biochemical methods to detect insufficient levels of cortisol followed by methods to establish the cause. Of such methods (which include the ACTH stimulation test, insulin induced hypoglycemia test and x-ray exams) the autoantibody test for the presence of antibodies to adrenal cortex is of prime significance as over 70% of patients with Addison's disease have autoimmune etiology. Indirect immunofluorescence (IF) on human or primate adrenal-cortex provide a simple, precise, and reliable method of detecting autoantibodies in Addison's disease<sup>3-6</sup>. Adrenal cortex antibodies (AcAb) are present in greater than 90% of patients with recent-onset autoimmune Addison's disease. AcAb are also markers of potential Addison's disease.

## PRINCIPLES OF PROCEDURES

Antibodies to the adrenal cortex are detected by indirect immunofluorescence (IF) using a primate adrenal substrate. In this IFA method, patient serum is incubated on tissue sections to allow binding of antibodies to the specific antigens of the substrate. Any immunoglobulins and other serum proteins not bound to the tissue sections are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of AcAb is determined by a specific apple-green reaction of the adrenal cortical cell cytoplasm. The titer (the reciprocal of the highest dilution giving a positive reaction) is determined by testing serial dilutions<sup>7</sup>.

## PRODUCT INFORMATION

### Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

### Materials Provided

ImmuGlo™ Anti-Adrenal cortex Antibody Test Kit **REF** 1171

Kits contain sufficient reagents to perform 48 determinations each.

8 x	<b>SORB</b> <b>SLD</b> <b>6</b>	<b>6 well Primate Adrenal Cortex Substrate Slides</b>
1 x 0.5 ml	<b>CONTROL</b> <b>+</b> <b>ADR</b> *	<b>Adrenal Antibody Positive Control.</b> Contains human serum.
1 x 0.5 ml	<b>CONTROL</b> <b>-</b> *	<b>Negative Control.</b> Contains human serum.
1 x 5 ml	<b>IgG-CONJ</b> <b>FITC</b> *	<b>Anti-human IgG FITC Conjugate.</b> <b>Protect from light.</b>
1 x 60 ml	<b>BUF</b> *	<b>Buffered Diluent.</b>
2 vials	<b>BUF</b> <b>WASH</b>	<b>Phosphate Buffered Saline (PBS).</b> Dissolve each vial to 1 liter.
1 x 5.0 ml	<b>MOUNTING</b> <b>MEDIUM</b> *	<b>Mounting Medium.</b> Do not freeze.
1 x 1.0 ml	<b>EVANS</b>	<b>Evan's Blue Counterstain.</b>
1 x 12	<b>COVER</b> <b>SLD</b>	<b>Coverslips</b>

\* Contains < 0.1% NaN<sub>3</sub>

### Materials Required But Not Provided

Fluorescence microscope  
Micropipette or Pasteur pipette  
Serological pipettes  
Staining dish (e.g. Coplin jar)  
Small test tubes (e.g. 13 x 75 mm) and test tube rack  
Distilled or deionized water  
1 liter container  
Wash bottle  
Paper towels  
Incubation chamber

### WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>14</sup>.

**WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

## SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### Test Method

#### A. Screening

1. Dilute each patient serum 1:4 with the Sample Diluent provided (0.1 ml serum + 0.3 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

## B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:4. The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

Number six tubes 1 through 6. Add 0.3 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

<b>Tubes</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
<b>Serum</b>	0.1 ml			
	+			
<b>Buffered Diluent</b>	0.3 ml	0.2 ml	0.2 ml	0.2 ml
<b>Transfer</b>		0.2 ml	0.2 ml	0.2 ml
<b>Final dilution</b>	1:4	1:8	1:16	1:32 etc.

## QUALITY CONTROL

Both a positive and negative control serum should be included with each test run. The negative control should show no specific fluorescence, whereas the positive control should have 2+ or greater staining intensity.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

## RESULTS

The anti-adrenal cortex antibody test should be reported as negative, or positive with titer. Read for specific staining of the adrenal cortical cell cytoplasm.

Various other tissue antibodies, such as anti-nuclear (ANA), anti-mitochondrial (AMA), anti-smooth muscle (ASMA), anti-endomysial antibodies (EMA) may also be observed on primate adrenal gland substrate. Sera exhibiting these reactions should be reported negative for anti-adrenal antibodies (see photo 1 at the end of this document).

## LIMITATIONS OF THE PROCEDURE

Different staining patterns can be observed. Most sera stain the whole cortex with brighter apple-green fluorescence of the glomerulosa zone. A few sera react only on the fasciculate and reticularis zones.

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
## REFERENCES • REFERENCIAS • LITERATUR • RIFERIMENTI

1. Muir A, Schatz DA, Maclaren NK. Autoimmune Addison's Disease. Springer semin Immunopathol 1993; 14:275-284.
2. Nerup J. Addison's Disease—clinical studies: a report of 108 cases. Acta Endocrinol 1974; 76:127
3. Betterle C, Volpato M, Rees Smith R et al. Adrenal cortex and steroid 21-hydroxylase autoantibodies in adult patients with organ-specific autoimmune diseases: Markers for low progression to clinical addison's disease. J Clin Endocrinol and Metab; 1997, 82:932-938.
4. Betterle C, Volpato M, Rees Smith R et al. Adrenal cortex and steroid 21-hydroxylase autoantibodies in children with organ-specific autoimmune diseases: Markers for high progression to clinical addison's disease. J Clin Endocrinol and Metab; 1997, 939-942
5. Laureti S, DeBellis A, Igino Muccitelli V et al. Levels of adrenocortical autoantibodies correlate with the degree of adrenal dysfunction in subjects with preclinical addison's disease. J Clin Endocrinol and Metab; 1997, 83:3507-3511.
6. Betterle C, Volpato M, Pedini B et al. Adrenal-cortex autoantibodies and steroid-producing cells autoantibodies in patients with addison's disease: comparison of immunofluorescence and immunoprecipitation assays. J Clin Endocrinol and Metab; 1999, 84:618-622.
7. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
8. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health. [HHS Pub. No. (CDC) 1993, 93-8395].

**Photo 1. Adrenal cortex positive for anti-adrenal antibodies.**

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**IMMCO**  
**DIAGNOSTICS**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)