# ImmuGlo™ Anti-Myelin Associated Antibody (anti-MAG) IFA

IVD

PRODUCT INSERT

REF

Code 1172

48 Determinations

## INTENDED USE

For the screening and detection by indirect immunofluorescence (IF) of primarily anti-myelin associated glycoproteins (anti-MAG) and other glycolipid autoantibodies in human serum.

## SUMMARY AND EXPLANATION

Autoimmune responses of the peripheral nervous system, recognized as *peripheral neuropathies*, are manifestations associated with autoantibodies against various neural glycoconjugates. These neuropathies can be acute, chronic, involve axonal degeneration, or demyelination. Autoimmune neuropathies can be further divided into monoclonal gammopathies and polyclonal inflammatory polyneuropathies like *Guillain-Barré syndrome*, *Chronic Inflammatory Demyelinating Polyneuropathy (CIDP)*, *multifocal motor neuropathy (MMN)*, and *paraneoplastic neuropathies*. In these diseases, there is a significant overlap of the involved auto-antigens which mediate the pathogenic mechanism. The following peripheral nerve specific autoantibodies are found in these neuropathies:

- a) **anti-Myelin associated glycoprotein (MAG),**
- b) **anti-acidic glycolipids like sulfoglucoronyl paragloboside (SGPG),**
- c) **anti-gangliosides**
- d) **anti-compact myelin associated proteins like P0, P2, and peripheral myelin protein 22 (PMP22)<sup>1-6</sup>.**

The epitope (HNK1) recognized by the human MAG autoantibodies is a sulfated oligosaccharide. This same epitope is shared by SGPG, P0 and PMP22<sup>7</sup>. Neuropathies associated with anti-MAG with IgM paraproteinemia are usually a heterogenous disease group, slowly progressive with evidence of demyelination and a variable degree of axonal loss usually associated with gait ataxia. Of all peripheral neuropathy cases with IgM paraproteinemia, 50% possess anti-MAG antibodies<sup>8</sup>. It is perceived that these autoantibodies might interfere with the process of myelination, with myelin maintenance, or with axon-Schwann cell interactions. Hence, the detection of these autoantibodies is useful for the clinician, as it suggests active demyelination in a peripheral neuropathy.

The IF technique is a sensitive method for the screening and detection of these anti-nerve myelin associated proteins and ganglioside autoantibodies. If the specimen yields no immunoreactivity by immunofluorescence the result is reported as negative. The specificity of the specimen, determined positive by IFA, should be confirmed by Western immunoblotting (anti-MAG) or by ELISA (anti-gangliosides).

## PRINCIPLES OF PROCEDURE

Antibodies to MAG and other peripheral nerve glycoconjugates are detected by indirect immunofluorescence using a primate peripheral nerve substrate. Patient sera are incubated on the section to allow binding of antibodies to the specific antigens of the substrate. Any immunoglobulins and other serum protein not bound to the tissue sections are removed by rinsing. Bound antibodies of the IgM class are detected by incubation of the substrate with fluorescein-labeled antihuman IgM conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of antibodies of a specificity is determined by their specific reactions to the nerve tissue. The anti-MAG antibodies

specifically stain the periaxonal oligodendroglial membrane, the Schwann cell periaxonal membranes, and the inner and the outer mesaxons of the nerve sheath.

## PRODUCT INFORMATION

### Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

### Materials Provided

ImmuGlo™ anti-MAG IFA **REF** 1172

Kits contain sufficient reagents to perform 48 determinations each.

8 x **SORB** **SLD** **6**

**6 well primate peripheral nerve substrate slides**

1 x 0.5 ml **CONTROL** + **MAG** \*

**MAG Positive Control**

1 x 0.5 ml **CONTROL** - \*

**Negative Control.** Contains human serum.

1 x 5 ml **IgM-CONJ** **FITC** \*

**Goat antihuman IgM FITC Conjugate.**  
Ready to use. **Protect from light.**

1 x 60 ml **BUF** \*

**Sample Diluent.**

2 vials **BUF WASH**

**Phosphate Buffered Saline (PBS).** Dissolve each vial to 1 liter.

1 x 5.0 ml **MOUNTING MEDIUM** \*

**Mounting Medium.** Do not freeze.

1 x 1.0 ml **EVANS**

**Evan's Blue Counterstain.**

1 x 12 **COVER SLD**

**Coverslips.** Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

### Materials Required But Not Provided

Fluorescence microscope  
Micropipette or Pasteur pipette  
Serological pipettes  
Staining dish (e.g. Coplin jar)  
Small test tubes (e.g. 13 x 75 mm) and test tube rack  
Distilled or deionized water  
1 liter container  
Wash bottle  
Paper towels  
Incubation chamber

### WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>22</sup>.

**WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

## SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### Test Method

#### A. Screening

1. Dilute each patient serum **1:10** with the Buffered Diluent provided (0.1 ml serum + 0.9 ml Diluent). **Do not** dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for **10-15 minutes**. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) of the **Negative Control** to well #1. Similarly apply **1 drop of Positive Control** to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur pipette, apply **1 drop** of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides **30 minutes** at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use a wash bottle. Transfer slide immediately into a Coplin jar filled with PBS and wash **10 minutes**. Repeat process with all remaining slides.
8. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the **Conjugate** dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) to each well. Repeat process with all remaining slides.
9. Replace the lid on the incubation chamber. Incubate **30 minutes** at room temperature.
10. Remove a slide from the incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for **10 minutes**. Repeat process with all remaining slides. If desired, **2-3 drops** of Evans blue counterstain may be added to the final wash.  
**NOTE:** Improper washing may lead to increased background fluorescence.
11. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **While slide is still wet mount the coverslip.** Place **3 drops** of the Mounting Medium evenly spaced on a coverslip placed on a paper towel and invert the slide onto the coverslip. To remove any air bubbles gently apply pressure along the edge of the coverslip. Avoid any movement of the coverslip. Repeat process with all remaining slides
12. Examine for specific fluorescence under a fluorescence microscope at a magnification of **200x** or greater.

Slides may be read as soon as prepared. However, because of the presence of an antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed. Slides should be stored in the dark at 2- 8°C.

**B: End Point Determination (Titration)**

A serum positive in the screening test may be further tested following **Steps 4 through 15** to determine the titer. Each test run should include the Positive and Negative Controls. Make serial twofold dilutions starting at 1:10 (see below). Using one slide, a serum may be tested at dilutions ranging from 1:10 to 1:320. If positive at a 1:320 dilution, the titer is reported as greater or equal to 320. Additional slides may be used to obtain endpoints for those sera still positive at a 1:320 dilution. The reciprocal of the highest dilution producing a positive reaction is the titer.

**Preparation of Serial Dilutions**

Number six tubes 1 through 6. Add 0.9 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6
Serum	0.1 ml					
	+					
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Transfer	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320 etc.

**Quality Control**

Both Positive and Negative Controls should be included with each test run. The negative control should show no specific fluorescence on the nerve tissue, whereas the positive control should have 1+ or greater staining intensity.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control.
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

**RESULTS**

Results of the test for anti-nerve antibodies should be reported as negative (<10) or positive with titer. For anti-MAG antibodies, read for staining of the periaxonal oligodendroglial membrane, the Schwann cell periaxonal membranes, and the inner and the outer mesaxons of the nerve sheath as shown in Figure 1 at the end of this document.

**LIMITATIONS OF THE PROCEDURE**

In some cases, sera positive for anti-MAG antibodies may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be retested at higher dilutions and, if positive, antibody titers determined.

The presence of two or more antibodies in a serum which are reactive with the substrate may cause an interference in their detection by immunofluorescence resulting in either, a failure to detect anti-MAG antibodies, or suppression of the titer. In addition, other distinct reaction patterns on the peripheral nerve sections may also be observed in sera positive for anti-SGPG, anti-P0, anti-MOG and anti-PMP22. Such reactions may not be easily distinguishable from anti-MAG due to the fact that they share the HNK1 epitope. It is recommended that all positive samples be confirmed for antibody specificity by western immunoblotting techniques.

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
## REFERENCES • REFERENCIAS • LITERATUR • RIFERIMENTI

1. Quarles RH, and Weiss MD. Autoantibodies associated with peripheral neuropathy. *Muscle & Nerve*; 22:800-822, 1999.
2. Griffin J. Antiglycolipid antibodies and peripheral neuropathies: links to pathogenesis. *Prog Brain Res*; 101: 313-323, 1994.
3. Quarles RH. Glycoproteins of the myelin sheaths. *J. Mol. Neurosci*; 8:1-12, 1997.
4. Kanda T, Yoshino H, Ariga T et al. Glycosphingolipid antigens in cultured bovine brain microvascular endothelial cells: Sulfoglucuronosyl paragloboside as a target of monoclonal IgM in demyelinating neuropathy. *J. Cell Biol*; 126:235-246, 1994.
5. Hammer JA, O'Shannessy DJ, and De LM et. al. Immunoreactivity of PMP-22, P0, and other 19 to 28kDa glycoproteins in peripheral nerve myelin of mammals and fish with HNK1 and related antibodies. *J. Neurosci Res*; 35:546-558, 1993.
6. Baba H, Daune GC, Ilyas AA et. al. anti-GM1 ganglioside antibodies with differing fine specificities in patients with multifocal motor neuropathy. *J. Neuroimmunol*; 25:143-150, 1989.
7. Field MC, Wing DR, Dwek RA et. al. Detection of multisulfated N-linked glycans in the L2/HNK1 carbohydrate epitope expressing neural adhesion molecule P0. *J. Neurochem*; 58:993-1000, 1992.
8. Latov N. Pathogenesis and therapy of neuropathies associated with monoclonal gammopathies. *Ann. Neurol*; 37:S32-S42, 1995.
9. Biosafety in Microbiological and Biomedical Laboratories. Center for Disease Control, National Institute for Health. (HHS Pub. No {CDC} 93-8395) 1993.
10. Chassande B, Leger JM, Younes-Chennoufi AB et al. Peripheral neuropathy associated with IgM monoclonal gammopathy: Correlations between M-Protein antibody activity and clinical/electrophysiological features in 40 cases. *Mus and Nerve*; 55-62, 1998.
11. Meucci M, Baldini L, Cappellari A et al. Anti-Myelin associated glycoprotein antibodies predict the development of neuropathy in asymptomatic patients with IgM monoclonal gammopathy. *Ann Neurol*; 46:119-122, 1999.

**Photo 1. Positive anti-MAG Reaction on Peripheral Nerve Substrate**

**Table 1. Incidence of anti-Nerve (anti-MAG) Antibodies in Neuropathies<sup>10</sup>**

<b>Disorder</b>	<b>n</b>	<b>positive</b>	<b>% positive</b>
Polyneuropathy+ IgM gammopathy	40	26	65
Demyelinating neuropathies	33	26	78
Axonal neuropathy	6	0	0

**Table 2. Association of anti-MAG Antibody Titers with Progression to Clinical Neuropathy<sup>11</sup>**

<b>anti-MAG Titer</b>	<b>no neuropathy n=17</b>	<b>subclinical n=7</b>	<b>confirmed clinical n=6</b>
<b>High 1:25,000,000 to 1:100,000</b>	1	3	3
<b>Low 1:6,400 to 1:200</b>	6	2	1
<b>Negative &lt;1:10</b>	10	2	2

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120**  
**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé  
EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
Tel (+31) 345 8570, Fax (+31) 346 7299  
[www.emergogroup.com](http://www.emergogroup.com)