# ImmuGlo™ Anti-Reticulin Antibody (ARA) Test System

For *in vitro* diagnostic use

**IVD**

CLIA Complexity: High  
CDC Analyte Identification Code: 0444  
CDC Test System Identification Code: 28331

## PRODUCT INSERT

Code 1115

48 Determinations

## INTENDED USE

An indirect immunofluorescence antibody test for the detection and semi-quantitation of anti-reticulin antibodies (ARA) in human serum.

## SUMMARY AND EXPLANATION

The detection of anti-reticulin antibodies (ARA) is performed primarily by indirect immunofluorescence. ARA was first described by Seah et al.<sup>1,2</sup> in 1971 in the sera of patients with gluten sensitive enteropathy: celiac disease (CD) or dermatitis herpetiformis (DH). ARA by definition are detected on rodent tissue. In 1973 Rizzetto and Doniach<sup>3</sup> identified five different ARA reaction patterns. Of these, the R1 pattern, characterized by peri-glomerular, peri-tubular, and peri-vascular staining of mouse or rat kidney is associated with gluten sensitive enteropathy. Both IgG and IgA immunoglobulin class ARA occur. The IgA class ARA is the more specific and sensitive marker of gluten sensitive enteropathy. The specificity of IgA-ARA has been reported to be 59-100% and its sensitivity between 30-95%<sup>4-11</sup>. IgG-ARA occur less frequently, either in conjunction with IgA class ARA or in CD patients who are IgA deficient<sup>11</sup>.

## PRINCIPLES OF PROCEDURES

In the indirect immunofluorescence method, patient serum is incubated on tissue sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgA and IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human immunoglobulin conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of ARA is demonstrated by an apple green fluorescence of the peri-tubular, peri-glomerular and peri-vascular fibers of the kidney. The titer (the reciprocal of the highest dilution giving a positive reaction) of the antibody is then determined by testing serial dilutions<sup>14</sup>

## PRODUCT INFORMATION

### Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

## Materials Provided

ImmuGlo™ ARA IFA REF 1115

Kits contain sufficient reagents to perform 48 determinations each.

8 x	<b>SORB</b> <b>SLD</b>	<b>6 well Rat Kidney Substrate Slides</b>
1 x 0.5 ml	<b>CONTROL +</b> <b>ARA</b> *	<b>ARA Positive Control.</b> Contains human serum.
1 x 0.5 ml	<b>CONTROL -</b> *	<b>Negative Control.</b> Contains human serum.
1 x 5 ml	<b>IgA/IgG-CONJ</b> <b>FITC</b> *	<b>Anti-human polyvalent FITC Conjugate.</b> Protect from light.
1 x 60 ml	<b>BUF</b> *	<b>Sample Diluent.</b>
2 vials	<b>BUF WASH</b>	<b>Phosphate Buffered Saline (PBS).</b> Dissolve each vial to 1 liter.
1 x 5.0 ml	<b>MOUNTING MEDIUM</b> *	<b>Mounting Medium.</b> Do not freeze.
1 x 1.0 ml	<b>EVANS</b>	<b>Evan's Blue Counterstain.</b>
1 x 12	<b>COVER</b> <b>SLD</b>	<b>Coverslips</b>

\* Contains < 0.1% NaN<sub>3</sub>

## Materials Required But Not Provided

Fluorescence microscope  
Micropipette or Pasteur pipette  
Serological pipettes  
Staining dish (e.g. Coplin jar)  
Small test tubes (e.g. 13 x 75 mm) and test tube rack  
Distilled or deionized water  
1 liter container  
Wash bottle  
Paper towels  
Incubation chamber

## WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>14</sup>.

**WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.**

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

## SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### Test Method

#### A. Screening

1. Dilute each patient serum 1:2.5 with the Sample Diluent provided (0.2 ml serum + 0.3 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

## B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:2.5. The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

Number four tubes 1 through 4. Add 0.3 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.2 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

<b>Tubes</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
<b>Serum</b>	<b>0.2 ml</b>			
	<b>+</b>			
<b>Buffered Diluent</b>	<b>0.3 ml</b>	<b>0.2 ml</b>	<b>0.2 ml</b>	<b>0.2 ml</b>
				
<b>Transfer</b>		<b>0.2 ml</b>	<b>0.2 ml</b>	<b>0.2 ml</b>
<b>Final dilution</b>	<b>1:2.5</b>	<b>1:5</b>	<b>1:10</b>	<b>1:20 etc.</b>

## QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the peri-nuclear, peri-glomerular and peri-vascular reticulin fibers, whereas the Positive Control should have 2+ or greater staining intensity of the tubules of these structures.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

## RESULTS

The results of the tests for ARA should be reported as negative (<2.5), positive greater or equal to 20, or preferably, positive with titer.

Read for specific peri-glomerular, peri-tubular and peri-vascular reticulin fiber staining. Other detectable antibodies include anti-smooth muscle antibodies (ASMA), antinuclear antibodies (ANA) and anti-mitochondrial antibodies (AMA).

## LIMITATIONS OF THE PROCEDURE

In some cases, sera positive for ARA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titer determined.

The presence of two or more antibodies in a serum which are reactive with the same tissue may cause an interference in their detection by immunofluorescence. This interference may cause either a failure to detect ARA or a suppression of its titer if the interfering antibody has a higher titer than ARA. The most common cause of the interference phenomenon in ARA tests

is the coexistence of ASMA. It is recommended that patients sera which also contain ASMA be tested further at higher dilutions.

In some patients with CD and DH, ARA tests may be negative. In such cases repeat testing may yield additional positive results. A negative test for ARA does not rule out a diagnosis of GSE. CD and DH patients on a gluten free diet may have low ARA titers or even be negative.

In some patients with celiac disease and IgA deficiency, the IgA-class ARA antibodies are absent. However, such patients are usually positive for IgG class ARA.

## EXPECTED VALUES

As seen in Table 1 at the end of this document, ARA are highly specific markers for GSE. The presence of ARA seems to be related to the intestinal pathology both in celiac disease and dermatitis herpetiformis rather than to the skin lesions in the latter.

## PERFORMANCE CHARACTERISTICS

The ImmuGlo™ Anti-Reticulin Antibody (ARA) Test kit was compared with another commercially available kit using a polyvalent conjugate and primate smooth muscle as a substrate. The comparison included a total of 100 sera: 50 from patients with clinically suspected celiac disease and 50 from normal controls. Sera were tested according to the procedure recommended by the manufacturer. The results were as follows:

		IMMCO ARA		
		POS (+)	NEG (-)	TOT (=)
<b>Other/Autre/Otro</b>	POS (+)	43	7	50
<b>Andere/Altro/Otra</b>	NEG (-)	0	50	50
<b>ARA IFA</b>	TOT (=)	43	57	100

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa:	93%
relative sensitivity/sensibilit� relative/sensibilidad relativa/relative Sensitivit�t/sensibilit� relativa/sensibilidade relativa:	86%
relative agreement/concordance/concordancia relativa/relative �bereinstimmung/concordanza relativa/acordo relativo:	100%

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8 x	<b>SORB</b>   <b>SLD</b>   <b>6</b>
1 x 0,5 ml	<b>CONTROL</b>   <b>+</b>   <b>ARA</b> *
1 x 0,5 ml	<b>CONTROL</b>   <b>-</b> *
1 x 5 ml	<b>IgA/IgG-CONJ</b>   <b>FITC</b> *
1 x 60 ml	<b>BUF</b> *
2 flaconcini	<b>BUF</b>   <b>WASH</b>
1 x 5,0 ml	<b>MOUNTING</b>   <b>MEDIUM</b> *
1 x 1,0 ml	<b>EVANS</b>
1 x 12	<b>COVER</b>   <b>SLD</b>

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

## REFERENCES • REFERENCIAS • LITERATUR • RIFERIMENTI

1. Seah PP, Fry L, Hoffbrand AV and Holborow EJ. Tissue antibodies in dermatitis herpetiformis and adult celiac disease. *Lancet*; 1971, 1:834-836.
2. Seah PP, Fry L, Holborow EJ, et. al. Anti-reticulin antibody. Incidence and diagnostic significance. *Gut*; 1973, 14:311-315.
3. Rizzetto M and Doniach D. Types of "reticulin" antibodies detected in human sera by immunofluorescence. *J Clin Path*, 1973, 26:841-851.
4. Khoshoo V, Bhan MK, Unsworth DJ, Kumar R and Walker-Smith JA. Anti-reticulin antibodies: useful adjunct to histopathology in diagnosing celiac disease, especially in a developing country. *J Pediatr Gastroenterol Nutr*; 1988, 7:864-866.
5. Kumar V, Hemedinger E, Chorzelski TP, Beutner EH, Valeski JE and Kowalewski C. Reticulin and endomysial antibodies in bullous diseases-comparison of specificity and sensitivity. *Arch Dermatol*; 1987, 123:1179-1182.
6. Lancaster-Smith M, Kumar P, Clark ML, Marks R and Johnson GD. Anti-reticulin antibodies in dermatitis herpetiformis and adult celiac disease. Their relationship to a gluten free diet and jejunal histology. *Br J Dermatol*; 1975, 92:37-42.
7. Ljunghall K, Scheynius A, Jonsson J, Schilling W and Forsum U. Gluten-free diet in patients with dermatitis herpetiformis. Effect on the occurrence of antibodies to reticulin and gluten. *Arch Dermatol*; 1983, 119:970-974.
8. Rossi TM, Kumar V, Lerner A, Heitlinger LA, Tucker N and Fisher J. Relationship of endomysial antibodies to jejunal mucosal pathology: specificity towards both symptomatic and asymptomatic celiacs. *J. Pediatr Gastroenterol Nutr*; 1988, 7:858-863.
9. Vainio E, Kosnai I, Hallstom O, Karpati S, Maki M and Reunala T. Antigliadin and antireticulin antibodies in children with dermatitis herpetiformis. *J Pediatr Gastroenterol Nutr*; 1986, 51:375-379.
10. Hallström O. Comparison of IgA-class reticulin and endomysium antibodies in celiac disease and dermatitis herpetiformis. *Gut*; 1989, 30:1225-1232.
11. Kumar V, Lerner A, Valeski JE, Beutner EM, Chorzelski TP and Kowalewski C. Relative sensitivity and specificity of serodiagnostic tests to gluten sensitive enteropathy. A comparative studies on IgA class endomysial (IgA-EMA) and IgA-class anti-reticulin antibodies (IgA-ARA) in celiac disease. In "Serologic Diagnosis of Celiac Disease", Chorzelski TP, Beutner EH, Kumar V and Zalewski TK, Eds, CRC Press, Boca Raton; 1990, 127-136.
12. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York; 1987, 3rd Ed, 3-40.
13. Biosafety in Microbiological and Biomedical Laboratories Centers for Disease Control, National Institutes of Health; 1993, (HHS Pub. No{CDC} 93-8395).

**Photo 1. Positive ARA staining reaction on rat kidney, 200X. Note staining of the reticulin fiber.**

**Photo 2. Negative ARA staining reaction.**

**Table 1. Incidence of ARA**

<b>Clinical Condition</b>	<b>No. Pos/Total</b>	<b>% Positive</b>
Confirmed Celiacs		
On Gluten	32/38	84
On GFD	4/37	11
Suspected Celiacs		
On Gluten	18/30	60
On GFD	2/30	7
Dermatitis Herpetiformis	19/54	35
Family Member of CD	2/76	3
Disease Controls		
Ulcerative Colitis	0/11	0
Crohn's Disease	0/74	0
Pemphigus	0/38	0
Pemphigoid	0/60	0
Blood Donors	0/113	0

GFD = Gluten Free Diet; adapted from Kumar et al. <sup>11</sup>

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive  
Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**IMMCO**  
**DIAGNOSTICS**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé  
EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,  
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)