


1.


2.


3.


4.


5.

6.

7.

A5341.4
March '04

MELISA™ TAB

Micropin enzyme-linked immunosorbent assay for thyroid antibodies

Combined Anti-Tg/TPO

Combined anti-thyroglobulin antibodies and anti-thyroid peroxidase antibodies assay kit

Code: M2496

Anti-Tg

Anti-thyroglobulin antibodies assay kit

Code: M2196

Anti-TPO

Anti-thyroid peroxidase antibodies assay kit

Code: M2396

Instructions for Use

Gebrauchsinformation

English: pages 1 to 13

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
Cambridge Life Sciences
Cambridgeshire Business Park, Angel Drove

Ely, Cambs. CB7 4DT UK


Tel: 01353 645200 Fax: 01353 645250

e-mail: clssales@byk.co.uk

www.cambridgelifesciences.co.uk

CAMBRIDGE
LIFE
SCIENCES

1. Intended Use

The Combined Anti-Tg/TPO, Anti-Tg and Anti-TPO MELISA™ TAB products are sandwich enzyme immunoassays for the quantitative detection of antibodies to thyroglobulin (anti-Tg) and to thyroid peroxidase (anti-TPO) in human serum.

The anti-Tg standards are calibrated against the WHO 1st Reference Preparation (1971) MRC 65/93.

The anti-TPO standards are calibrated against the reference preparation NIBSC 66/387.

2. Background

Autoimmune thyroid gland disorders are characterised by the detection of anti-thyroid antibodies against Tg and TPO antigens. TPO has been identified as the specific antigenic determinant of the thyroid microsomal antigen.

The presence of auto-antibodies to thyroid antigens correlates with the degree of lymphocytic infiltration of the thyroid gland, the hallmark of Hashimoto's thyroiditis. Circulating auto-antibodies to Tg and TPO are present in the serum of over 90% of patients with thyroid autoimmune diseases such as Hashimoto's thyroiditis, Graves' disease, idiopathic hypothyroidism and sub-acute thyroiditis. In each case, the prevalence of anti-TPO antibodies is greater than that of anti-Tg antibodies.

It is reported that low levels of autoimmune antibodies predict at-risk pregnancy. Furthermore, a report using an EIA test for anti-thyroglobulin and anti-recombinant TPO demonstrated a 100% increase in the rate of spontaneous miscarriage in women who had detectable serum thyroid auto-antibodies in their first trimester of pregnancy.

Thyroid auto-antibodies are detected using immunoassays such as passive haemagglutination, indirect fluorescence antibody (IFA), enzyme immunoassay (EIA) and radioimmunoassay (RIA) techniques. MELISA™ TAB Anti-Tg uses human thyroglobulin. MELISA™ TAB Anti-TPO uses recombinant human thyroid peroxidase which does not contain contaminating thyroglobulin and/or mitochondria found in other microsomal antigen preparations. Both of these assays are in an EIA test format.

3. Principle

MELISA™ TAB products employ a unique antigen-coated micropin technology, which is ideal for the batch-screening of large numbers of samples for anti-thyroid antibodies. The methods utilise a non-competitive sandwich enzyme immunoassay system.

First Incubation

MELISA™ TAB micropins are provided coated with purified antigen (human thyroglobulin and recombinant thyroid peroxidase respectively). When the micropins are immersed into wells containing standards, controls or diluted sera, any anti-Tg or anti-TPO antibodies present will bind to the micropin surface. The micropins are then removed and washed in wash buffer.

15. Literatur

1. Salvi, M. *et al* (1988) Role of auto-antibodies in the pathogenesis and association of endocrine autoimmune disorders; *Endocrine Reviews* **9** (4): 450-465
2. Volpe, R. (1977) The role of autoimmunity in hypoendocrine and hyperendocrine function. *Ann Int Med* **87**: 86-99
3. Yoshida, H. *et al* (1978) Association of serum antithyroid antibodies with lymphocytic infiltration of the thyroid gland: studies of seventy autopsied cases. *J Clin Endo Metab* **46** (6): 859-862
4. Cowchock, S. *et al* (1984) Subclinical autoimmune disease and unexplained abortion. *Am. J Obstet Gynecol* **150**: 367-371
5. Maier, D.B. *et al* (1989) Subclinical autoimmunity in recurrent aborters. *Fertil Steril* **51** (2): 280-285
6. Stagnaro-Green, A. *et al* (1990) Detection of at-risk pregnancy by means of highly sensitive assays for thyroid auto-antibodies. *JAMA* **264** (11): 1422-1425
7. Centers for Disease Control/National Institutes of Health Manual: Biosafety in Microbiological and Biomedical Laboratories (1984)
8. Hijmans, W. *et al* (1961) Serological overlap between lupus erythematosus, rheumatoid arthritis and thyroid auto-immune disease. *Brit Med J* **2**: 909-914
9. Roman, S.H. *et al* (1986) Thyroid auto-antibodies in HLA-genotyped type I diabetic families: sex-limited DR5 association with thyroid microsomal antibody. *Clin Endocrin* **25**: 23-33
10. Beever, K. *et al* (1989) Highly sensitive assays of auto-antibodies to thyroglobulin and to thyroid peroxidase *Clin Chem* **35** (9): 1949-1954
11. Stites, D.P. *et al* (1994) Basic and clinical immunology. 8th edition. Appleton & Lange.

Second Incubation

The micropins are then dipped into wells containing anti-IgG-horseradish peroxidase (HRP) conjugate, which will bind to any captured antibodies. Unbound conjugate is removed by washing in wash buffer.

Third Incubation

The micropins are then dipped into wells containing a colourless substrate. The reaction is stopped with the supplied stop solution. The intensity of the yellow colour formed is proportional to the concentration of anti-thyroid antibody bound in the first incubation.

4. Kit Contents

5 standards* (1.0mL ready to use)

Standard	Anti-Tg IU/mL	Anti-TPO IU/mL
1	0	0
2	300	100
3	750	250
4	3000	1000
5	7500	2500

(*For Combined Anti-Tg/TPO MELISA™ TAB, product code:M2496, 2 x 5 standards are included i.e. 1 set Anti-Tg and 1 set Anti-TPO).

- 1 vial wash buffer concentrate (x25), 20mL
- 1 vial sample diluent, 100mL
- 1 vial conjugate (anti-IgG-HRP), 12mL
- 1 vial substrate (TMB*), 21mL
- 1 vial stop solution (oxalic acid), 12mL
- 1 foil sachet containing 1 set of 96 antigen-coated micropins and 1 flat-well plate
- 1 vial Positive control (1.0mL ready to use)
- 1 vial Negative control (1.0mL ready to use)
- 2 round-well plates
- 2 wash trays
- 1 instruction leaflet
- 1 QC certificate

* TMB = 3,3', 5,5'-tetramethyl benzidine.

5. Storage



The kit should be stored at 2-8°C. Do not use the reagents beyond their expiry date. Do not freeze. Keep all reagents away from direct sunlight.

6. Safety Precautions

For *in-vitro* diagnostic use only.

For Professional use only.

The standards and controls contain human source material. Although found negative when tested for HIV-1 and HIV-2 antibodies, HCV and hepatitis B surface antigen, no test can guarantee their absence. Therefore, the standards and controls should be handled using the same safety precautions employed when handling any potentially infectious material.

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Used standards, samples and controls, pipette tips and plates should be handled as clinical waste and incinerated or disposed of in accordance with local rules. Other reagents should be diluted and flushed down the drain. It is recommended that gloves are worn when handling such items.

Safety data sheets are available on request.

7. Sample Handling

The assay may be performed on human serum samples. Samples should be assayed within 24 hours of collection or stored frozen at -15°C or colder. Repeated freeze-thawing is not advisable. Do not heat treat samples.

8. Additional Reagents and Equipment

Deionised water.
Precision micropipettes and disposable tips to deliver 5 - 1000µL.
Multichannel micropipette or repeating dispenser to deliver 100 and 200µL.
500mL measuring cylinder for reagent preparation.
96-well microplate reader with 450nm filter.
Software package (optional).

9. Procedural Precautions

Carefully read instructions for use before starting with the assay.

Allow all reagents to equilibrate to room temperature (18-25°C) before use for a minimum of two hours.

Avoid the use of icteric, lipaemic or grossly haemolysed samples.

Always change pipette tips between different standards, samples or controls to prevent sample carryover.

Never allow the same pipette tip to be used with different reagents. Special care is needed to prevent contamination of the substrate by the conjugate.

The substrate should be colourless. Any coloration indicates substrate contamination and the substrate should be discarded.

The micropin washing procedure is critical for the successful performance of the test, especially between conjugate and substrate incubations (i.e. the second and third incubations). Do not use the kit beyond the expiry date given on the box.

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
10. Assay Procedure

For Combined Anti-Tg/TPO MELISA™ TAB plate layout, see diagram below. Product Code: M2496.

In the Combined Anti-Tg/TPO assay kit, anti-Tg is measured on the left hand side of the plate and the anti-TPO on the right hand side.

1. Prepare the wash buffer as follows: dilute contents of the **wash buffer concentrate** (x25) vial to 500mL with deionised water.
2. Dilute patient samples 1/200 using the **sample diluent** e.g. 5µL sample added to 995µL diluent. The **standards** and **controls** do not require dilution.
3. Prepare the sample plate by dispensing 100µL of each **standard**, diluted patient sample or **controls** into appropriate wells of a round-well plate. It is recommended that samples be tested in duplicate.
4. Prepare for the micropipette washing stage by pouring wash buffer into each of the two wash trays, so that it fills the central reservoir.

5. Remove the **antigen-coated micropins** from the foil sachet. Transfer the micropins to the sample plate and incubate for 10 minutes at room temperature (18-25°C). During the incubation period, prepare the conjugate plate by dispensing 100µL **conjugate** into each well of the second round-well plate.
6. After incubation, transfer the micropins to the first wash tray and gently agitate by sideways and up and down movements for at least 10 seconds. Flick off the excess wash buffer or gently drain on absorbent paper. Repeat in the second wash tray, afterwards ensuring that the micropins are free of wash buffer and bubbles.
7. Transfer the micropins to the conjugate plate and incubate for 10 minutes at room temperature (18-25°C). During the incubation period, prepare the substrate plate by dispensing 200µL of **substrate** into each well of the flat-well plate. Also, prepare for the final micropin washing stage by replenishing the wash trays with wash buffer so that it fills the central reservoir.
8. After incubation, wash the micropins as in step 6.
9. Transfer the micropins to the substrate plate and incubate for 10 minutes at room temperature (18-25°C).
10. After incubation, carefully remove the micropins from the plate and discard them. Add 100µL of **stop solution** to each well. The developed colour will change from blue to yellow.
11. Measure the absorbance at 450nm on a 96-well microtitre plate reader.

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11. Calculation of Results

For each assay, prepare a calibration curve by plotting mean absorbance against standard concentration on linear graph paper, and interpolate unknowns. Alternatively, use a computerised curve-fit program. Any sample giving values above the standard range should be further diluted and re-tested.

12. Quality Control

Good laboratory practice requires that quality control samples be included in every run to check on assay performance. High and low controls for anti-Tg and anti-TPO antibodies are provided within the kits. These controls are ready-to-use and do not require dilution. If either control value falls outside the quoted range on the QC Certificate, the results are invalid and the assay should be repeated.

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
2. Minimum detectable concentration

The minimum detectable concentration, defined as the concentration equal to 2 standard deviations from the mean of the sample diluent, was found to be:



17.5 IU/mL (anti-Tg)
1.8 IU/mL (anti-TPO)

3. Reference Ranges

MELISA™ TAB was used to determine the anti-Tg/TPO levels of 100 serum samples measured in duplicate from normal blood donors with no apparent abnormalities. The data was evaluated and the following ranges obtained:

	Anti-Tg (IU/mL)	Anti-TPO (IU/mL)
Normal range	<164.3	<72.5
Borderline	164.3 - 217.6	72.5 - 97.6
Positive	>217.6	>97.6

It is advised that each laboratory establishes its own reference range.

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