EN

IMMCO
DIAGNOSTICS

Anti-Mitochondrial M2 Antibody (AMA) ELISA

IVD

PRODUCT INSERT

REF 1163 Anti-Mitochondrial (M2) Antibody ELISA 96 Determinations

INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of anti-mitochondrial antibodies (AMA) in human serum. The presence of AMA can be used as an adjunct to clinical and other laboratory findings in the diagnosis of primary biliary cirrhosis.

SUMMARY AND EXPLANATION

Primary biliary cirrhosis (PBC) is a disease of the liver characterized by inflammatory obliteration of the intrahepatic bile ducts^{1,2}. PBC is characterized by the presence of anti-mitochondrial antibodies (AMA). AMA detected by indirect immunofluorescence (IF) occur in approximately 90% of patients with PBC and in 3-10% of patients with autoimmune hepatitis as well as in certain other disorders such as syphilis and drug induced hepatitis³. As many as nine different staining reactions have been described by indirect IF, but only one, the M2 reaction is specific for PBC⁴. Inability to easily distinguish different AMA reactions compromises the utility of indirect IF test methods for diagnosing PBC.

Biochemical and immunochemical studies have identified PBC specific mitochondrial antigen to be localized on the inner mitochondrial membrane and is identified to be the mitochondrial enzyme pyruvate dehydrogenase (PDH). Subsequent studies identified PBC antigen to be 70 kD molecular weight. The mitochondrial antigen associated with PBC is the E2 subunit of a functionally related enzyme family of 2-oxo-acid dehydrogenase complex (2-OADC)⁵⁻⁷.

The Immco™ anti-mitochondrial M2 antibody ELISA is a specific and sensitive immunoassay for the detection of AMA associated with PBC.

PRINCIPLES OF PROCEDURES

The test is performed as a solid phase enzyme labeled immunosorbent assay (ELISA) in microwells coated with purified mitochondrial antigen. Controls, calibrators and patient serum samples are incubated in the microwells allowing anti-mitochondrial antibodies present in the serum to bind the antigen. Unbound antibody and other serum proteins are removed by washing the microwells. Bound antibodies are incubated with an enzyme labeled anti-human IgG conjugate. Unbound conjugate is removed by washing the microwells. Specific enzyme substrate (pNPP) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of the substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in Enzyme Units per milliliter (EU/ml).

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.**

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20°-25°C) prior to use.

When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered

EN

potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials⁸.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use beyond expiration date on the label.

Materials provided

ImmULisa™ Mitochondrial (M2) Antibody ELISA **REF** 1163

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	MICROPLATE M2	Microplate with individual breakaway microwells coated with mitochondrial antigen
1 x 1.5 ml	CALIBRATOR A M2 *	Ready to use Calibrator A (<i>green cap</i>). Human serum containing mitochondrial antibodies.
1 x 1.5 ml	CALIBRATOR B M2 *	Ready to use Calibrator B (<i>violet cap</i>). Human serum containing mitochondrial antibodies.
1 x 1.5 ml	CALIBRATOR C M2 *	Ready to use Calibrator C (<i>blue cap</i>). Human serum containing mitochondrial antibodies.
1 x 1.5 ml	CALIBRATOR D M2 *	Ready to use Calibrator D (<i>yellow cap</i>). Human serum containing mitochondrial antibodies.
1 x 1.5 ml	CONTROL + M2 *	Ready to use Positive Control (<i>red cap</i>). Contains human serum positive for mitochondrial antibodies.
1 x 1.5 ml	CONTROL - *	Ready to use Negative Control (<i>whitecap</i>). Contains human serum.
1 x 12 ml	IgG-CONJ ALKPHOS *	Ready to use anti-human Alk. Phos. Conjugate . Color coded pink.
1 x 60 ml	DIL *	Ready to use Serum Diluent . Color coded blue.
1 x 12 ml	SUBSTRATE *	Ready to use Enzyme Substrate . Contains pNPP. Protect from light.
1 x 12 ml	STOP	Ready to use Stop Solution .
2 x	BUF WASH	Powder Wash Buffer . Reconstitute to one liter each.

* Contains <0.1% NaN₃

Symbols used on labels:

LOT Lot number

REF Catalog number

EN

-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer or Automatic Microplate washer capable of dispensing 200 µl
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes


- Before starting with the assay read carefully the product insert.
- Let serum specimens and test reagents equilibrate at room temperature before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Good washing technique is critical. If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.

EN


Test Method

- Step 1** Let all reagents and specimens equilibrate at room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use only the Ready to Use Low Calibrator D (*vial with yellow cap*).
or For a **semi-quantitative determination** use the Ready to Use Calibrators A through D as depicted in the sample layout below.

QUALITATIVE DETERMINATION

SEMI-QUANTITATIVE DETERMINATION

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **0.5 ml** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder .
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples to the appropriate microwells as per protocol sheet.
Note: Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within 1 hour from adding Stop Solution.
- Step 15** Read absorbance of each microwell at **405 nm** using a single or 405/630nm dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <20 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

RESULTS

Calculations


The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

2. SEMI-QUANTITATIVE DETERMINATION

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such serum sample should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor.

EN

Interpretation

The following serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values. These may vary with the population examined.

AMA value	Interpretation
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

LIMITATIONS OF THE PROCEDURE

The results obtained serve only as an aid in the diagnosis of PBC and should not be interpreted as diagnostic in themselves. Serum of some PBC patients may be negative for AMA.

EXPECTED VALUES

Incidence of AMA by ELISA

Disease	Incidence %
Primary Biliary Cirrhosis	95-100
Autoimmune Hepatitis	0
Systemic Lupus Erythematosus	0
Rheumatoid Arthritis	0
Systemic Scleroderma	0

Note: The incidence of AMA may vary depending upon the patient population studied and the criteria of selection.

PERFORMANCE CHARACTERISTICS

The Immulisa™ anti-mitochondrial antibody test was compared with a commercially available indirect immunofluorescence test kit for the detection of anti-mitochondrial antibodies in human serum.

A total of 113 sera were obtained from a clinical reference laboratory. They were identified as positive or negative for AMA by indirect immunofluorescence and were tested according to the procedures recommended by the manufacturer. The results are summarized in the following table:

Comparison of ELISA and IFA for detecting AMA				
		Positive	Negative	Total
Other ELISA	Positive	35	2	37
	Negative	2	74	76
	Total	37	76	113

Agreement: 96.5 %
Sensitivity: 94.6 %
Specificity: 97.4 %

EN

Precision:

Based on 10 replicates, the intra-assay and inter-assay coefficient of variation (CV) were calculated and found to be between 5.4 and 9.5%, depending upon the concentration of the antibody:

	mean value (EU/ml)	inter-assay %CV	intra-assay %CV
Sample 1	146	5 %	7 %
Sample 2	82	10 %	7 %

Recovery:

Samples with known AMA concentrations were mixed with appropriate dilutions of another positive sample with known amounts. AMA levels of the mixed samples were determined and from the values obtained the percent recovery calculated:

	Ab. conc added (EU/ml)	Ab. conc. obtained (EU/ml)	% Recovery
Sample 1	106	108	102
Sample 2	88	78	89
Sample 3	53	56	105

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Kaplan MM. Primary biliary cirrhosis, Review. *N Engl J Med*; 1987, 316:521-528.
2. Colucci G, Schaffner F and Paronetto F. In situ characterization of the cell-surface antigens of the mononuclear cell infiltrates and bile duct epithelium in primary biliary cirrhosis. *Clin Immunol & Immunopathol*; 1986, 41:35-.
3. Patrick SC, Leung, Ross L et al. Antimitochondrial antibodies in primary biliary cirrhosis. *Seminars Liv Dis*; 1997, 17:61-69.
4. Berg PA, Klein R, Lindenborn-Fotinos J et al. Heterogeneity of anti-mitochondrial antibodies. *Lancet*; 1982, 2:1423-1426.
5. Yeaman SJ, Fussey SP, Danner DJ, and Bassendine MF. Primary biliary cirrhosis: identification of two major M2 mitochondrial autoantigens. *Lancet*; 1988, 2:1067-1070.
6. Van de Water J, Gershwin ME, Leung P et al. The autoepitope of the 74-kD mitochondrial autoantigen of primary biliary cirrhosis corresponds to the functional site of dihydrolipoamide acetyltransferase. *J Exp Med*; 1988, 167:1791-1799.
7. Fussey SM, Ali ST, Guest JR et al. Reactivity of primary biliary cirrhosis sera with *Escherichia coli* dihydrolipoamide acetyltransferase (E2p); characterization of the main immunogenic region. *Proc Natl Acad Sci USA*; 1990, 87:3987-3991.
8. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
9. N. Zurgil, R. Bakimer, M. Kaplan et al. Anti-pyruvate dehydrogenase autoantibodies in primary biliary cirrhosis. *J Clin Immunol*; 1991, 5:239-245.

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com