

## COPPER Multi-Purpose (MPR) Liquid Reagent

### KIT SPECIFICATIONS:

Cat. No.	Quantity	Reagent	Storage
GL553CU	2 x 60 ml	COPPER	04 – 22 °C
	1 x 10 ml	COPPER - Standard	
GL563CU	6 x 60 ml	COPPER	04 – 22 °C
	1 x 10 ml	COPPER - Standard	

### INTENDED USE:

In Vitro Diagnostic reagent pack for the determination of Copper in serum and plasma on automated and semi-automated analysers.

### PRINCIPLE OF THE TEST: <sup>1</sup>

Copper forms with 4-(3, 5-dibromo-2-pyridylazo)-N-ethyl-N-sulfopropylanine, a chelate complex. The increase of absorbance of this complex can be measured and is proportional to the concentration of total copper in the sample.

### WARNINGS AND PRECAUTIONS:

*For In Vitro Diagnostics Use Only - For Professional Use Only.*

Carefully read and follow instructions for use. Deviations from the described procedure may alter performance of the assay.

### Components Colour and Appearance:

Reagent 1: Light pink liquid.

Refer to Laboratory's QC program for actions to be taken. In case of serious damage to the bottle and/or cap, resulting in product leakage and/or contamination, do not use the reagent pack and contact your distributor.

### Safety Precautions:

**CAUTION:** Take all necessary precautions required when handling laboratory reagents. Material Safety Data Sheet is available upon request.

### Handling precautions:

- Do not use components past the expiry date stated on the Bottles.
- Do not Freeze Reagents.
- Do not use components for any purpose other than described in the "Intended Use" section.
- Do not interchange caps among components as contamination may occur and compromise test results.
- Refer to local legal requirements for safe waste disposal.

### INSTRUMENTS:

Instrument application procedures are available upon request.

### COMPONENT COMPOSITION:

Component	Ingredients	Concentration in Tests
Reagent 1	Acetate Buffer pH 5.0	0.2 mol/l
	4-(3,5-dibromo-2-pyridylazo) -N-ethyl-N-sulfopropylaniline	0.02 mmol/l
Standard	Copper	200 µg/dl
		(31.46 µmol/l)

### REAGENT PREPARATION AND STABILITY:

**Reagent 1 and standard** are ready to use.

Before use, mix reagent by gently inverting each bottle.

If stored and handled properly, components are stable until expiry date stated on label.

### TYPE OF SPECIMEN: 2

Use serum or heparin plasma as specimen.

It is recommended to follow NCCLS procedures (or similar standardised conditions) regarding specimen handling. Specimen should be collected in an appropriate sample container, with proper specimen identification. Serum/plasma should be separated from cells within 2 hours after collection.

*Stability:* up to 48 hours at 2-8°C.

### TEST PROCEDURE:

Materials required but not supplied:

Description	Catalogue No.
General Chemistry Control Level 1	GL922
General Chemistry Control Level 2	GL932
General Laboratory Equipment	N/A
Saline solution 0.9 g/l NaCl	N/A
Photometer	N/A

### Assay procedure:

Wavelength: 580 nm

Temperature: 37°C

Optical path: 1 cm light path

	Blank	Calibrator	Sample
<b>Reagent 1</b>	1 ml	1 ml	1 ml
<b>Sample</b>	----	----	50 µl
<b>Calibrator</b>	----	50 µl	----
Gently mix and Incubate at 37°C for 5 minutes. Read Optical Density (OD).			

### Calibration:

Using recommended calibrator or standard provided, calibrate the assay:

- When using a new reagent kit or changing lot number.
- Following preventive maintenance or replacement of a critical part of the photometer used.
- When Quality Control results are out of range.

### Quality Control:

All clinical laboratories should establish an Internal Quality Control program. Verify instrument and reagent performance with recommended controls or similar. The values obtained for QC should fall within manufacturer's acceptable ranges or should be established according to the Laboratory's QC program.

Controls should be assayed:

- Prior reporting patient results.
- Following any maintenance procedure on the photometer used.
- At intervals established by the laboratory QC Programme.

### CALCULATION:

$$\text{Concentration of Copper} = \frac{OD_{\text{Sample}}}{OD_{\text{Calibrator}}} \times \text{Concentration of Calibrator}$$

(Conversion factor: Qty in µmol/l = Qty in µg/dl x 0.157)

### Expected Values: 1

In serum	Adult Man	70 – 140 µg/dl	11.0 – 22.0 µmol/l
	Adult Women	80 – 155 µg/dl	12.6 – 24.3 µmol/l

Each laboratory should establish its own reference range. Copper results should always be reviewed with the patient's medical examination and history.

### PERFORMANCE CHARACTERISTICS:

Performance evaluation can vary with the instrument used. Data obtained in each individual laboratory may differ from results obtained in-house.

### Linearity:

This assay is linear up to 500 µg/dl (78.5 µmol/l).

For samples with a higher concentration, dilute 1:1 with 0.9% NaCl (9g/l) and re-assay. Multiply result by 2.

### Precision:


















Test 15	Mean (mmol/l)	Max (mmol/l)	Min (mmol/l)	SD	% CV
Level 1	32.0	33.2	30.9	0.61	1.90

### BIBLIOGRAPHY:

- Abe A., Yamashita S., Noma A., Clin Chem., 552-554, 35 (1989)
- C. A. Burtis, E.R. Ashwood. Tietz Fund. Of Clin. Chem. 5<sup>th</sup> ed. 30:54 and 973

### SYMBOLS:

The following symbols are used in the labelling of Glenbio systems:

	In Vitro Diagnostics		Catalogue No
	Batch Code		Content
	Reagent		Antibody
	Calibrator		Substrate
	Buffer		Aqueous Standard
	CE Mark - Device complies with the Directives 98/79/EC		
	Storage temperature		Reconstitute with
	Expiry Date (Last day of the month)		Manufactured By
	Biological risk		Consult Instruction for Use

Manufactured By: GLENBIO LTD,  
10 Kilbegs Road, Antrim, Co. Antrim, BT41 4NN, UK



Tel/Fax: +44(0)2879659842  
Email: info@glenbio.com  
Web: www.glenbio.com