Fr. Juin 2019

REF HM 059-2X1

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REF HM 059-2X1

INTENDED USE

This reagent, prepared from human serum (heathy donor), is intended for use as an unassayed precision control with **DOT, ELISA, FIDIS and IFA** procedures for determination of autoantibodies.

DESCRIPTION

IMMUNO-TROL reagents have been designed for use with **DOT, ELISA, IFA and FIDIS** methods to monitor the precision of laboratory testing procedures for the analytes of autoimmunity

This **IMMUNO-TROL Neg** reagent lot is designed to produce a negative reaction for the research of autoantibodies.

Each laboratory should establish its own target value for each marker and for each new lot on a minimum of 30 assays.

PACKAGING

2X1 : 2 vials of 1 mL

Green Caps

CONTROL -

Results must be inferior to the cut off of the technique used for the search of the parameters.

WARNING

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatite B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

- This reagent contain <0.06% ProClin®300 and <0.1% of sodium azide. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

- **IMMUNO-TROL NEG reagents have been developed according to EC regulation n°1272/2008.**

-  At this concentration, ProClin®300 is a skin sensitizer, prolonged or repeated exposure may cause allergic reaction in certain sensitive individuals.

H317: May cause an allergic skin reaction

PROCEDURE

IMMUNO-TROL reagents should be analyzed in the same manner as patient samples, according to the instructions supplied by manufacturers of the kits used.

STABILITY AND STORAGE

- Store the reagents at +2°C to +8°C in their original package.
- Do not use reagents beyond the expiration date.
- Open vial stability : 60 days at +2°C to +8°C.
- Do not freeze.

LIMITS

- Performances and levels of reactivity of the **IMMUNO-TROL** reagents may vary with different manufacturer's tests and assay procedures.

- **IMMUNO-TROL** reagents must not be substituted for the positive and the negative control reagents provided with manufacturer's tests.

- Because **IMMUNO-TROL** reagents are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.


- **IMMUNO-TROL** reagents are provided for quality assurance and must not be used for calibration.

- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IMMUNO-TROL** reagents for its particular use and to establish guidelines for interpretation of results obtained with **IMMUNO-TROL** reagents.

- Adverse storage conditions or use of outdated reagents may produce erroneous results.

- **IMMUNO-TROL NEG** reagents have been tested and showed a negative reactivity in test procedure of the following  kits: DOT, ELISA, IFA and FIDIS.

SYMBOLS USED

	Biological risk		Temperature limitation
	Read instructions for use	REF	Catalog Number
	In Vitro Diagnostic Use	LOT	Lot Number
	EC Declaration of Conformity		Use by
	Manufacturer	CONTROL -	Negative control