Fr. Juillet 2019

**REF** HM 064-2X05

**REF** HM 064-5X1

**REF** HM 064-5X05

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**REF** HM 064-5X05

**INTENDED USE**

This reagent, prepared from human serum, is intended for use as an unassayed precision control with **DOT, ELISA and FIDIS** procedures for determination of the following autoantibodies:

|   | Codes           |
|---|-----------------|
| • human and animal <b>anti-Rheumatoid Factors</b> IgM isotype | RF (H&A)<br>IgM |
| • <b>anti-cyclic citrullinated peptide</b> IgG isotype        | CCP IgG         |

**DESCRIPTION**

**IMMUNO-TROL** reagents have been designed for use with **DOT, ELISA and FIDIS** methods to monitor the precision of laboratory testing procedures for the analytes listed in this insert.

Each **IMMUNO-TROL** reagent lot is designed to produce a positive reaction for each autoantibodies above listed.

**IMMUNO-TROL** reagents do not have assigned values.

Each laboratory should establish its own target value for each marker and for each new lot on a minimum of 30 assays.

**PACKAGING**

|                         |                      |
|-------------------------|----------------------|
| 2X05 : 2 vials of 0,5mL | 5X1 : 5 vials of 1mL |
| 5X05 : 5 vials of 0,5mL | <b>CONTROL +</b>     |

**Blue caps**

The antigenic specificities are confirmed on the vial label. The antibody codes are specified above.

**WARNING**

- Human sources for the preparation of this control have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus or hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.
- The reagent contains <0.06% of ProClin®300 and <0.1% of sodium azide. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.
-  At this concentration, ProClin®300 is a skin sensitizer, prolonged or repeated exposure may cause allergic reaction in certain sensitive individuals.

The **IMMUNO-TROL** reagents have been developed according to EC regulation n°1272/2008.

H317: May cause an allergic skin reaction

The **IMMUNO-TROL** reagents have been developed according to EC regulation n°1272/2008.

H317: May cause an allergic skin reaction

**PROCEDURE**

**IMMUNO-TROL** reagents should be analyzed in the same manner as patient samples, according to the instructions supplied by manufacturers of the test kits used. If the value of one parameter is out of the kit standard range, the reagent may be diluted following the kit's instruction for use.


**STABILITY AND STORAGE**

- Store the reagents at +2°C to +8°C when not in use.
- Do not use reagents beyond the expiration date.
- Open vial stability : 60 days at +2°C to +8°C.
- Do not freeze.

**LIMITS**

- Performances and levels of reactivity of the **IMMUNO-TROL** reagents may vary with different manufacturer's test and assay procedures.
- **IMMUNO-TROL** reagents must not be substituted for the positive and the negative control reagents provided with manufacturer's tests.
- Because **IMMUNO-TROL** reagents are formulated with human sera and non human derived constituents, incompatibilities between this reagent and test kits may exist.
- **IMMUNO-TROL** reagents are provided for quality assurance and must not be used for calibration.
- It is the responsibility of each laboratory to implement its own quality assurance program, to determine the suitability of **IMMUNO-TROL** reagents for its particular use and to establish guidelines for interpretation of results obtained with **IMMUNO-TROL** reagents.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- **IMMUNO-TROL** reagents have been tested and showed a reactivity in test procedure of the following  kits: DOT, ELISA and FIDIS.

**SYMBOLS USED**

|   |                              |   |                        |
|---|------------------------------|---|------------------------|
|  | Biological risk              |  | Temperature limitation |
|  | Read instructions for use    | <b>REF</b>  | Catalog Number         |
|  | In Vitro Diagnostic Use      | <b>LOT</b>  | Lot Number             |
|  | EC Declaration of Conformity |  | Use by                 |
|  | Manufacturer                 | <b>CONTROL +</b>  | Positive control       |
| Changes made on the previous version  |                              |   |                        |