

PARACETAMOL

KEY POINTS:


- **Paracetamol (Acetaminophen) 2 Reagent Assay Kit**
- **Easy Enzyme reconstitution**
- **Adaptable to multiple clinical chemistry platforms**
- **No interference from up to 1g/L N-acetyl cysteine or Paracetamol metabolites**
- **Serum and plasma samples**

INDICATION:


Quantitative measurement of paracetamol as an aid to monitor and diagnose paracetamol overdose.

PRINCIPLE OF THE ASSAY:

The method is based on the use of an enzyme specific for the amide bond of acylated aromatic amines. It cleaves the paracetamol molecule, yielding p-aminophenol, which reacts specifically with o-cresol in ammoniacal copper solution to produce a blue colour. The assay is specific for the parent compound and does not detect paracetamol metabolites.



ASSAY SCHEME:





REFERENCE RANGE:

Paracetamol result	Result Interpretation
66 - 199 µmol/L (10 - 30 mg/L)	Therapeutic level
>1990 µmol/L (>300 mg/L)	Toxic at 4 hours after ingestion
>330 µmol/L (> 50 mg/L)	Toxic at 12 hours after ingestion

ORDERING INFORMATION

PARACETAMOL (Acetaminophen)
ASSAY KIT

  98/79/EC
For in vitro diagnostic use and
Professional Use only.

CONFIGURATION

5 x Lyophilised Enzyme (**R1**)
5 x 20mL Enzyme Diluent (**R1**)
1 x 100mL Colour Reagent (**R2**)
2 x 0.5mL 0mg/L Sera Calibrator
2 x 0.5mL 302mg/L Sera Calibrator
1 x 3mL 302mg/L Aqueous Calibrator
(option)

PART NUMBER


K8003

Cambridge Life Sciences Ltd.
14 St. Thomas' Place, Cambridgeshire Business Park, Ely, Cambs, CB7 4EX, UK
Telephone: +44 (0)1353 645200
E-mail: sales@clsdiagnostics.com
<http://www.clsdiagnostics.com/>



PERFORMANCE:

Accuracy



Linearity

0.02 - 5.00 mmol/L (3 - 756 mg/L) using sera calibrators

K8003 Paracetamol method performance on Pentra C400 was compared to the external quality scheme (NEQAS) results

Precision

Sample		Level 1	Level 2	Level 3
N		46	46	46
Mean (mmol/L)		0.0857	0.2869	0.9003
Mean (mg/L)		12.94	43.32	135.95
Repeatability	SD	0.0023	0.0032	0.0075
	%CV	2.68	1.12	0.83
Between Day	SD	0.0007	0.0039	0.0090
	%CV	0.82	1.36	1.00
Between Run	SD	0.0027	0.0040	0.0091
	%CV	3.15	1.39	1.01
Total	SD	0.0032	0.0060	0.0137
	%CV	3.73	2.09	1.52


Interference

Interferent	Concentration
Ascorbic Acid	1.76g/L (10.0mmol/L)
Total Bilirubin (unconjugated)	300.0mg/L (513µmol/L)
Direct Bilirubin (conjugated)	300.0mg/L (513µmol/L)
Haemoglobin	5.0g/L (77.6µmol/L)
Triglycerides	10.0 g/L (11.3mmol/L)
Intralipids	1% (w/v)
N-acetyl cysteine	1 g/L

Recovery

Target (mmol/L)	Mean (mmol/L)	Difference (mmol/L)	Target (mg/L)	Mean (mg/L)	Difference (mg/L)	Recovery %
0.31	0.313	0.003	46.87	47.32	0.44	100.9
0.62	0.632	0.012	93.74	95.49	1.75	101.9
1.26	1.252	-0.008	190.51	189.03	-1.48	99.2
2.52	2.522	0.002	381.02	380.89	-0.13	100.0
3.00	3.052	0.052	453.60	458.81	5.21	101.6
5.00	5.152	0.028	756.00	779.00	23.00	103.0

On-board Reagent Stability



Calibrator Standardisation

Paracetamol calibrators are manufactured using primary calibration material, Acetaminophen (99.5% - 100.5%) that meets USP specifications. They are manufactured gravimetrically and tested against independent controls.