
# LISA TRACKER

## Duo Certolizumab

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### DEFINITION

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**Arends S & al.** The formation of autoantibodies and antibodies to TNF- $\alpha$  blocking agents in relation to clinical response in patients with ankylosing spondylitis. *Clin Exp Rheumatol.* 2010 Sep-Oct;28(5):661-8.

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**Bendtzen K.** Is There a Need for Immunopharmacologic Guidance of Anti-Tumor Necrosis Factor Therapies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 867–870.

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**Desroches M & al.** Treatment failure with antagonists of TNF- $\alpha$ : mechanisms and implications for the care of patients. *Eur. Cytokine Netw.*, Vol. 21 n° 4, December 2010, 226-31.

**DeVries & al.** Inefficacy of infliximab in ankylosing spondylitis is correlated with antibody formation. *Ann Rheum Dis.* 2007 Jan;66(1):133-4.

**Edrees & al.** Anti-tumor necrosis factor (TNF) therapy in rheumatoid arthritis: correlation of TNF-alpha serumlevel with clinical response and benefit from changing dose or frequency of infliximab infusions. *Clin Exp Rheumatol.* 2005 Jul-Aug;23(4):469-74.

**Jaminitski & al.** The presence or absence of antibodies to infliximab or adalimumab determines the outcome of switching to etanercept. *Rheum Dis.* 2011 Feb;70(2):284-8.

**Koren & al.** Recommendation on risk-based strategies for detection and characterization of antibodies against biotechnology products. *Journal of Immunological Methods*, 333 (2008) 1-9.

**Korswagen LA & al.** Venous and Arterial Thromboembolic Events in Adalimumab-Treated Patients With Anti-adalimumab Antibodies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 877–883.

**Krzysiek R & al.** Circulating Concentration of Infliximab and Response to Treatment in Ankylosing Spondylitis: Results From a Randomized Control Study. *Arthritis & Rheumatism (Arthritis Care & Research)* Vol. 61, No. 5, May 15, 2009, pp 569–57.

**Marotte & al.** Circulating tumour necrosis factor-alpha bioactivity in rheumatoid arthritis patients treated with infliximab: link to clinical response. *Arthritis Res Ther.* 2005;7(1):R149-55.

**Mire-Sluis & al.** Recommendation for the design and optimization of immunoassays used in the detection of host antibodies against biotechnology products. *Journal of Immunological Methods*, 289 (2004) 1-16.

**Piketny ML & al.** Biotin: an emerging analytical interference, *Ann Biol Clin* 2017 ; 75(4) : 366-8.

**Radstake TRDJ & al.** Formation of antibodies against infliximab and adalimumab strongly correlates with functional drug levels and clinical responses in rheumatoid arthritis. *Ann Rheum Dis* 2009;68:1739–1745.

**Ram Raj Singh, & al.** TNF $\alpha$  blockade in human diseases : mechanisms and future. *Clin.Immunol.*, 2008 February ; 126(2):121-136.

**Seow CH & al.** Trough serum infliximab: a predictive factor of clinical outcome for infliximab treatment in acute ulcerative colitis. *Gut.* 2010 Jan;59(1):49-54.

**Sparado A & al.** Switching from infliximab or etanercept to adalimumab in resistant or intolerant patients with spondyloarthritis: a 4-year study. *Rheumatology (Oxford).* 2010 Jun;49(6):1107-11.

**Steenholdt C, & al.** Cut-off levels and diagnostic accuracy of infliximab trough levels and anti-infliximab antibodies in Crohn's disease. *Scandinavian Journal of Gastroenterology*, March 2011, Vol. 46, N°3, Pages 310-318.

**Takeuchi & al.** Baseline tumour necrosis factor alpha levels predict the necessity for dose escalation of infliximab therapy in patients with rheumatoid arthritis. *Ann Rheum Dis.* 2011 Jul;70(7):1208-15.

**Ternant D & al.** An enzyme-linked immunosorbent assay for therapeutic drug monitoring of Infliximab. *Ther Drug Monit.*, 2006 April ; 28(2):169-174.

**Van den Bemt BJJ & al.** Anti-infliximab antibodies are already detectable in most patients with rheumatoid arthritis halfway through an infusion cycle: an open-label pharmacokinetic cohort study. *BMC Musculoskeletal Disorders* 2011.

**Vetterlein & al.** No antibodies to PEG detected in patients treated with certolizumab pegol. *Ann Rheum Dis* 2008;67(Suppl II):236.

**Wolbink GJ & al.** Development of Anti-infliximab Antibodies and Relationship to Clinical Response in Patients With Rheumatoid Arthritis. *ARTHRITIS & RHEUMATISM* Vol. 54, No. 3, March 2006, pp711–715.

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
**REF** LTC 005

**42 determinations**
**DEFINITION**

**LISA-TRACKER Duo Certolizumab** () is an enzyme linked immunoassay (ELISA) for the quantitative determination of Certolizumab pegol (anti-TNF $\alpha$ ) and anti-Certolizumab pegol antibodies in human serum samples. These tests can be separately or simultaneously done by following the standardized assay protocols.

**DIAGNOSTIC VALUE**

Anti-TNF $\alpha$  are therapeutic agents widely used to treat patients with various inflammatory diseases. Certolizumab pegol is one of the anti-TNF $\alpha$  recommended for the treatment of the rheumatoid arthritis. This Fab' part of humanized monoclonal antibody conjugated to polyethylene glycol is able to bind TNF $\alpha$ . It blocks the action of TNF $\alpha$  responsible for the inflammatory state. However, during the treatment, some patients can develop antibodies against Certolizumab pegol. Consequently, the plasmatic level of anti-TNF $\alpha$  decreases and simultaneously the disease symptoms reappear or increase.

**LISA-TRACKER Duo Certolizumab** () allows the detection of 2 parameters: Certolizumab pegol and anti-Certolizumab pegol antibodies. This kit allows the physician to monitor the level of these 2 parameters in patient sera.

**SAMPLES COLLECTION AND HANDLING**

- The test should be performed on plasma (EDTA, heparin or citrate collection tubes) or on serum.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- To avoid any non-specific binding, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

**ASSAY PRINCIPLE**
**A. Dosage of Certolizumab pegol**

The TNF $\alpha$  is coated onto a polystyrene microtiter plate (6 strips of 8 wells).

- First, the diluted sample is added to the TNF $\alpha$  coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Biotinylated anti-polyethylene glycol antibodies (anti-PEG-biot) are added. After incubation, unbound antibodies are removed by washing
- Then horseradish peroxidase labelled streptavidin is added. The streptavidin binds to the complex formed with biotinylated anti-PEG antibodies. After incubation, the wells are washed again to eliminate any excess of conjugate.

- The bound enzyme is revealed by addition of substrate TMB (3,3',5,5' tetramethylbenzidin). The colour intensity is proportional to the amount of Certolizumab pegol.
- Adding H<sub>2</sub>SO<sub>4</sub> allows to stop the enzymatic reaction.
- After stopping the reaction by H<sub>2</sub>SO<sub>4</sub>, the optical density is read by a spectrophotometer at 450nm.

A range of calibration allows to define the quantity of Certolizumab pegol of each patient samples expressed in  $\mu$ g/mL.

**B. Dosage of anti-Certolizumab pegol**

The Certolizumab pegol is coated onto a polystyrene microtiter plate (6 strips of 8 wells).

- First, the diluted sample is added to the antibody coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Biotinylated Certolizumab pegol is added. After incubation, unbound antibodies are removed by washing
- Then horseradish peroxidase labelled streptavidin is added. The streptavidin binds to the complex formed with biotinylated Certolizumab pegol. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate TMB (3,3',5,5' tetramethylbenzidin). The colour intensity is proportional to the amount of anti-Certolizumab pegol antibodies.
- Adding H<sub>2</sub>SO<sub>4</sub> allows to stop the enzymatic reaction.
- After stopping the reaction by H<sub>2</sub>SO<sub>4</sub>, the optical density is read by a spectrophotometer at 450nm.

A range of calibration allows to define the quantity of anti-Certolizumab pegol antibodies of each patient samples expressed in AU/mL.

**REAGENTS**

3 reagent families :

Color	Certolizumab pegol reagents	anti-Certolizumab pegol antibodies reagents	Common reagents
cap of vials	Blue	White or colorless	Green, White, Black or Purple
microwells	Blue	colorless	-

### A) Specific reagents for the Certolizumab pegol determination

Strips of individual breakaway blue wells coated with human TNF $\alpha$ <b>MP</b>	6 strips
5 vials of « Certolizumab pegol » Standards, ( $\mu\text{g/mL}$ ). <b>CER CAL n</b> <u>Ready to use.</u> The vials can be reused several times. <b>The quantity of Certolizumab pegol is indicated on the vial label</b> <b>Blue caps</b>	5 x 1.5mL
«Positive control - Certolizumab pegol», ( $\mu\text{g/mL}$ ). <u>To dilute.</u> The vials can be reused several times. <b>CER CONTROL +</b> <b>The quantity of Certolizumab pegol is indicated on the vial label</b> <b>Blue cap</b>	1 x 250 $\mu\text{L}$
Biotinylated antibody vial. <u>Ready to use</u> <b>CER Ab BIOT</b> <b>Blue cap</b>	1 x 7.5mL

### B) Specific reagents for the anti-Certolizumab pegol antibodies determination

Strips of individual breakaway colorless wells coated with Certolizumab pegol. <b>MP</b>	6 strips
5 vials of « anti-Certolizumab pegol » Standards, (AU/mL). <b>A-CER CAL n</b> <u>Ready to use.</u> The vials can be reused several times. <b>The quantity of anti-Certolizumab pegol is indicated on the vial label</b> <b>Colorless cap.</b>	5 x 1.5mL
« Positive control – anti-Certolizumab pegol », (AU/mL). <b>A-CER CONTROL +</b> <u>To dilute.</u> The vial can be reused several times. <b>The quantity of anti-Certolizumab pegol is indicated on the vial label.</b> <b>Colorless cap.</b>	1 x 1mL
Biotinylated antibody vial. <u>Ready to use.</u> <b>A-CER Ab BIOT</b> <b>White cap.</b>	1 x 7.5mL

### C) Common reagents

HRP labelled Streptavidin. <u>Ready to use.</u> <b>CONJ HRP</b> <b>Green cap.</b>	1 x 12mL
Phosphate-Tween Buffer pH 7,2 (10x) – <u>To reconstitute with distilled water.</u> <b>White cap.</b> <b>BUF WASH 10x</b>	1 x 100mL
Phosphate Buffer pH 7,2 (10x) – <u>To reconstitute with distilled water.</u> <b>White cap.</b> <b>BUF PHOS 10x</b>	1 X 100ml
Substrate (TMB). <u>Ready to use</u> <b>SUBS TMB</b> <b>Black cap.</b>	1 x 12mL
Stop solution - H <sub>2</sub> SO <sub>4</sub> (0.25 N). <u>Ready to use.</u> <b>SOLN STOP</b> <b>Purple cap</b>	1 x 15mL

### MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- microplate spectrophotometer with 450 nm filter
- 8 channel pipettes

### STABILITY AND STORAGE

- Store reagents and micro-wells at +2°C/+8°C in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag (closed securely) with the desiccant.
- Store all components immediately after use again at +2°C/+8°C.

### SETUP

Except the TDL, which can be prepare in advance, all reagents must be prepared extemporaneously.

#### 1. Dilution and Wash buffer (TDL)

- Dilute concentrated Phosphate-Tween Buffer 1/10 in distilled water.

**BUF WASH 10x**

- Shelf life : 3 months at +2°C/+8°C (avoid to use if signs of contamination or other visible changes occur).

*NB. If there are crystals in the concentrated solution, warm the bottle up to +37°C for 15 minutes before use.*

#### 2. Phosphate buffer

- Dilute concentrated Phosphate Buffer 1/10 in distilled water.

**BUF PHOS 10x**

- Shelf life : 3 month at +2°C/+8°C (avoid to use if signs of contamination or other visible changes occur).

*NB. If there are crystals in the concentrated solution, warm the bottle up to +37°C for 15 minutes before use.*

#### 3. Preparation of samples and positive controls

##### a. Samples

##### - Certolizumab pegol determination

- Dilute to 1/7000 in TDL

Ex: perform a 1/100 dilution followed by a 1/70 dilution (10 $\mu\text{L}$  sample + 990 $\mu\text{L}$  TDL, vortex, add 10 $\mu\text{L}$  of the 1/100 dilution into 690 $\mu\text{L}$  of TDL).

- Vortex vigorously.

##### - Anti-Certolizumab pegol determination

- Dilute to 1/2 in TDL

Ex : 130 $\mu\text{L}$  sample + 130 $\mu\text{L}$  TDL

- Vortex vigorously.

##### b. Positive controls

##### - Certolizumab pegol determination

- Dilute to 1/7000 in TDL

Ex: perform a 1/100 dilution followed by a 1/70 dilution (10 $\mu\text{L}$  sample + 990 $\mu\text{L}$  TDL, vortex, add 10 $\mu\text{L}$  of the 1/100 dilution into 690 $\mu\text{L}$  of TDL).

- Vortex vigorously.

##### - anti-Certolizumab pegol determination

- Dilute to 1/2 in TDL

Ex : 130 $\mu\text{L}$  positive control + 130 $\mu\text{L}$  TDL

- Vortex vigorously.

#### 4. Use of ready-to-use biotinylated antibodies.

- Estimate the amount required for handling and transfer to a tube.

#### 5. Use of ready-to-use HRP Streptavidin conjugate.

- Estimate the amount required for handling and transfer to a tube.

#### 6. Use of ready-to-use substrate.

- Estimate the amount required for handling and transfer to a dark tube.

### PRECAUTIONS

---

Unpack all reagents in order to let them warm at room temperature (+18°C/+25°C) at least half an hour before starting the test.

⚠ The temperature of the reagents can impact the final result.

Check that all plates are well drained after each wash.

Avoid to use reagents if signs of contamination or other visible changes occur.

*Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.*

Reagents in solution (except for substrate buffer and stop solution) contain <0.1% of sodium azide and <0.6% of ProClin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

 At this concentration, ProClin® 300 is irritating to eyes and skin, and may be detrimental if enough quantity is ingested. It is a skin sensitizer; prolonged or repeated exposure may cause allergic reaction in certain sensitive individuals.

**LISA-TRACKER Duo Certolizumab** () has been developed according to EC regulation n°1272/2008 relating to the classification, labelling and packaging of substances and mixtures.

**LISA-TRACKER Duo Certolizumab** () has been optimized for the use as describe in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

### METHOD

---

#### 1. Preparing the test

Use the work sheet to note the sample locations.

Set out:

- 5 "standard" wells
- 1 well for positive control
- 1 well for each sample

For a simultaneous testing of the 2 parameters, detach the exact number of wells needed. Return unused wells to plastic pouch provided in the kit, with the desiccant bag.

#### Remark :

If a dispensing/diluting device is used, place the specific wells in the following order : dosage of Certolizumab and dosage of anti-Certolizumab.

#### 2. Samples, positive controls and standards incubation

Add 100 µl of standards, diluted controls or samples.

Incubate for 60 minutes at room temperature (+18°C/+25°C).

*Wash step:*

Remove the content of the wells by rapid inversion.

**Wash 3 times with 300µL of TDL buffer.**

**Wash 3 times with 300µL of phosphate buffer.**

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

#### 3. Incubation of biotinylated antibodies

Add 100µL of specific biotinylated antibodies in identified wells.

Incubate for 60 minutes at room temperature (+18°C/+25°C).

*Wash step:*

Remove the content of the wells by rapid inversion.

Wash 3 times with 300µL of dilution and washing buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

#### 4. Incubation of Conjugate

Add 100µL of conjugate.

Incubate for 30 minutes at room temperature (+18°C/+25°C).

*Wash step:*

Remove the content of the wells by rapid inversion.

Wash 3 times with 300µL of dilution and washing buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

#### 5. Incubation of Substrate

Add 100µL substrate into each well.

Incubate for 15 minutes at room temperature (+18°C/+25°C). Protect from light.

#### 6. Stop of the reaction

Add 100µl of H<sub>2</sub>SO<sub>4</sub> to each well.

#### 7. Reading

Read the optical density of each well at 450nm within 30 minutes after stopping reaction.

### RESULTS AND INTERPRETATION

---

#### A. Dosage of Certolizumab pegol

- The OD of the standard 1 should be at least 0.8.
- The positive control value should be comprised into the range indicated on the vial label.
- Trace the standard curve (polynomial curve), plotting the units of the 5 standard points (µg/mL) along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).
- The Certolizumab pegol value can be directly read on the curve.
- Samples with values greater than that of standard 1 may be diluted to obtain a more precise result. The number of units should be multiplied by the selected dilution.

#### B. Dosage of anti-Certolizumab pegol

- The OD of the standard 1 should be at least 0.8.
- The positive control value should be comprised into the range indicated on the vial label.

- Trace the standard curve (polynomial curve), plotting the units of the 5 standard points (AU/mL) along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).
- The anti-Certolizumab pegol value can be directly read on the curve.
- Samples with values greater than that of standard 1 may be diluted to obtain a more precise result. The number of units should be multiplied by the selected dilution.

## CHARACTERISTICS AND PERFORMANCE OF THE TEST

### Limits of detection / threshold values

Estimated on healthy patient samples:

- 152 samples for Certolizumab pegol
- 142 samples for anti-Certolizumab pegol antibodies

Limit of detection <b>Certolizumab pegol</b> Calculated with the range standard of 0 – 84µg/mL	Limit of detection <b>Anti-Certolizumab pegol</b>
3 µg/mL >98 <sup>th</sup> percentile	5 AU/mL >99 <sup>th</sup> percentile

### Assay range

<b>Certolizumab pegol</b> Calculated with the range standard of 0 – 84µg/mL	<b>Anti-Certolizumab pegol</b>
3 µg/mL – 84 µg/mL	5 AU/mL - 160 AU/mL

### Interfering Substances Study


**LISA-TRACKER Duo Certolizumab** (Theradag) was evaluated to assess potential cross reactivity to other antibodies and interference from serum components (cryoglobulins, rheumatoid factors, heterophilic antibodies, high levels of triglycerides, bilirubin, IgG and/or IgM, and C1q proteins, autoantibodies). No interference was detected.

### Precision

Parameters	Intra-run *4 tests in a same assay **3 tests in a same assay		Inter-runs *4 tests 4 different assays **3 tests 6 different assays	
	Mean	CV (%)	Mean	CV (%)
<b>Certolizumab (µg/mL)**</b>	5.8	10.4	5.5	11.1
	28.8	16.5	28.1	12.9
	71.5	6.1	56.2	22.6
<b>Anti-Certolizumab (AU/mL)*</b>	10	6.7	10	9.5
	40	3.6	40	11.5
	85	4.4	84	12.6

### Correlation

The results of the Certolizumab (CZP) determination obtained with the extended standard range (0-84µg/mL) were compared to the results obtained with the previous standard range (0-12µg/mL). A linear regression analysis of the two results showed that both are equivalent.

## LIMITS

The presence of biotin in patients' specimens can potentially impact the immunodosages using the Streptavidin-Biotin technology.

## QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. **IMMUNO-TROL Certolizumab** or **anti-Certolizumab** controls (Theradag), Cat.n: LTC 002-PC or LTC 003-PC, contain **Certolizumab** or antibodies directed against **Certolizumab**. These materials are to be assayed in the same manner as the unknown samples.

## REFERENCES

**Ainsworth MA & al.** Tumor necrosis factor-alpha binding capacity and anti-infliximab antibodies measured by fluidphase radioimmunoassays as predictors of clinical efficacy of infliximab in Crohn's disease. *Am J Gastroenterol.* 2008 Apr;103(4):944-8.

**Afif W & al.** Clinical Utility of Measuring Infliximab and Human Anti-Chimeric Antibody Concentrations in Patients With Inflammatory Bowel Disease. *Am J Gastroenterol* 2010; 105:1133–1139.

**Arends S & al.** The formation of autoantibodies and antibodies to TNF-α blocking agents in relation to clinical response in patients with ankylosing spondylitis. *Clin Exp Rheumatol.* 2010 Sep-Oct;28(5):661-8.

**Baert F & al.** Influence of Immunogenicity on the Long-Term Efficacy of Infliximab in Crohn's Disease. *N Engl J Med* 2003;348:601-8.

**Bartelds GM & al.** Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA.* 2011 Apr 13;305(14):1460-8.

**Bentzen K.** Is There a Need for Immunopharmacologic Guidance of Anti-Tumor Necrosis Factor Therapies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 867–870.

**Candon S & al.** Clinical and biological consequences of immunization to infliximab in pediatric Crohn's disease. *Clin Immunol.* 2006 Jan;118(1):11-9.

**Choy & al.** Efficacy of a novel PEGylated humanized anti-TNF (CDP870) in patients with rheumatoid arthritis : a phase II double-blinded, randomized, dose escalating trial. *Rheumatology* 2002;41:1133-1137.

**Desroches M & al.** Treatment failure with antagonists of TNF- $\alpha$ : mechanisms and implications for the care of patients. Eur. Cytokine Netw., Vol. 21 n° 4, December 2010, 226-31.

**DeVries & al.** Inefficacy of infliximab in ankylosing spondylitis is correlated with antibody formation. Ann Rheum Dis. 2007 Jan;66(1):133-4.

**Edrees & al.** Anti-tumor necrosis factor (TNF) therapy in rheumatoid arthritis: correlation of TNF-alpha serumlevel with clinical response and benefit from changing dose or frequency of infliximab infusions. Clin Exp Rheumatol. 2005 Jul-Aug;23(4):469-74.

**Jaminitzki & al.** The presence or absence of antibodies to infliximab or adalimumab determines the outcome of switching to etanercept. Rheum Dis. 2011 Feb;70(2):284-8.

**Koren & al.** Recommendation on risk-based strategies for detection and characterization of antibodies against biotechnology products. Journal of Immunological Methods, 333 (2008) 1-9.

**Korswagen LA & al.** Venous and Arterial Thromboembolic Events in Adalimumab-Treated Patients With Anti-adalimumab Antibodies. ARTHRITIS & RHEUMATISM Vol. 63, No. 4, April 2011, pp 877–883.

**Krzysiek R & al.** Circulating Concentration of Infliximab and Response to Treatment in Ankylosing Spondylitis: Results From a Randomized Control Study. Arthritis & Rheumatism (Arthritis Care & Research) Vol. 61, No. 5, May 15, 2009, pp 569–57.

**Marotte & al.** Circulating tumour necrosis factor-alpha bioactivity in rheumatoid arthritis patients treated with infliximab: link to clinical response. Arthritis Res Ther. 2005;7(1):R149-55.

**Mire-Sluis & al.** Recommendation for the design and optimization of immunoassays used in the detection of host antibodies against biotechnology products. Journal of Immunological Methods, 289 (2004) 1-16.

**Piketny ML & al.** Biotin: an emerging analytical interference, Ann Biol Clin2017 ; 75(4) : 366-8.

**Radstake TRDJ & al.** Formation of antibodies against infliximab and adalimumab strongly correlates with functional drug levels and clinical responses in rheumatoid arthritis. Ann Rheum Dis 2009;68:1739–1745.

**Ram Raj Singh, & al.** TNF $\alpha$  blockade in human diseases : mechanisms and future. Clin.Immunol., 2008 February ; 126(2):121-136.

**Seow CH & al.** Trough serum infliximab: a predictive factor of clinical outcome for infliximab treatment in acute ulcerative colitis. Gut. 2010 Jan;59(1):49-54.

**Sparado A & al.** Switching from infliximab or etanercept to adalimumab in resistant or intolerant patients with spondyloarthritis: a 4-year study. Rheumatology (Oxford). 2010 Jun;49(6):1107-11.

**Steenholdt C, & al.** Cut-off levels and diagnostic accuracy of infliximab trough levels and anti-infliximab antibodies in Crohn's disease. Scandinavian Journal of Gastroenterology, March 2011, Vol. 46, N°3, Pages 310-318.

**Takeuchi & al.** Baseline tumour necrosis factor alpha levels predict the necessity for dose escalation of infliximab therapy in patients with rheumatoid arthritis. Ann Rheum Dis. 2011 Jul;70(7):1208-15.

**Ternant D & al.** An enzyme-linked immunosorbent assay for therapeutic drug monitoring of Infliximab. Ther Drug Monit., 2006 April ; 28(2):169-174.

**Van den Bemt BJF & al.** Anti-infliximab antibodies are already detectable in most patients with rheumatoid arthritis halfway through an infusion cycle: an open-label pharmacokinetic cohort study. BMC Musculoskeletal Disorders 2011.

**Vetterlein & al.** No antibodies to PEG detected in patients treated with certolizumab pegol. Ann Rheum Dis 2008;67(Suppl II):236.

**Wolbink GJ & al.** Development of Anti-infliximab Antibodies and Relationship to Clinical Response in Patients With Rheumatoid Arthritis. ARTHRITIS & RHEUMATISM Vol. 54, No. 3, March 2006, pp711–715.

## SUMMARY OF METHOD

### A) Sample Dilution

Certolizumab pegol	Anti-Certolizumab pegol
1/7000	1/2

### B) Positive Control Dilution

Certolizumab pegol	Anti-Certolizumab pegol
1/7000	1/2

### C) Procedure

Reagents	Procedure
Standards	100 $\mu$ L / well
Diluted positive controls	
Diluted samples	
Incubation	1 h at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300 $\mu$ L / well
<b>Washing*</b>	<b>Wash 3 times with phosphate buffer : 3 x 300<math>\mu</math>L / well</b>
Biotinylated antibodies	100 $\mu$ L / wells ( <b>specific reagents</b> )
Incubation	1 h at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300 $\mu$ L / well
HRP-Streptavidin	100 $\mu$ L / well
Incubation	30 minutes at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300 $\mu$ L / well
Substrate (TMB)	100 $\mu$ L / well
Incubation	15 minutes at room temperature. Protect from light.
Stop solution	100 $\mu$ L / well

\* Dry the microplate by tapping it gently on a towel to eliminate the excess of liquid.

**D) Configuration of the assays**

**a. 42 sera for Certolizumab and anti-Certolizumab pegol**

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Certolizumab  
pegol assay

anti-  
Certolizumab  
pegol assay

**b. 2 sera for Certolizumab and anti-Certolizumab pegol**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	Standard 5	Standard 5										
<b>B</b>	Standard 4	Standard 4										
<b>C</b>	Standard 3	Standard 3										
<b>D</b>	Standard 2	Standard 2										
<b>E</b>	Standard 1	Standard 1										
<b>F</b>	C+	C+										
<b>G</b>	Sera 1	Sera 1										
<b>H</b>	Sera 2	Sera 2										

Certolizumab  
pegol assay

anti-  
Certolizumab  
pegol assay

### c. 42 sera for anti-Certolizumab pegol determination

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