
LISA TRACKER

Duo Vedolizumab

REF LTV 005

Français

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DEFINITION

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BIBLIOGRAPHIE

Cherry LN & al., *Vedolizumab: an $\alpha 4\beta 7$ integrin antagonist for ulcerative colitis and Crohn's disease*, Therapeutic Advances in Chronic Disease, Vol.6(5), 224-233, 2015.

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Petrovic A. & al., *LPAM ($\alpha 4\beta 7$ integrin) is an important homing integrin on alloreactive T cells in the development of intestinal graft-versus-host disease*, Blood, 103, 4, 2004.

Piketty ML & al. Biotin: an emerging analytical interference, Ann Biol Clin 2017 ; 75(4) : 366-8.

Rosario M. and al., *Population pharmacokinetics-pharmacodynamics of vedolizumab in patients with ulcerative colitis and Crohn's disease*, Alimentary Pharmacology and Therapeutics, 2015.

Sandborn WJ. And al., *Vedolizumab as induction and maintenance therapy for crohn's disease*, The new England journal of medicine, 369, 711-721, 2013.

Soler D. and al., *The binding specificity and selective antagonism of vedolizumab, an anti- $\alpha 4\beta 7$ integrin therapeutic antibody in development for inflammatory bowel diseases*, The journal of pharmacology and experimental therapeutics, 330, 864-875, 2009.

Yu Jin and al., *Meta-analysis of the effectiveness and safety of vedolizumab for ulcerative colitis*, World Journal of Gastroenterology, 21, issue 20, 2015.

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
42 determinations
English

DEFINITION

LISA-TRACKER Duo Vedolizumab () is an enzyme linked immunoassay (ELISA) for the quantitative determination of Vedolizumab (anti- $\alpha 4\beta 7$ integrin) and anti-Vedolizumab antibodies in human serum samples. These quantitative tests can be separately or simultaneously done by following the standardized assay protocol.

DIAGNOSTIC VALUE

Vedolizumab is a humanized monoclonal antibody against $\alpha 4\beta 7$ integrin. It's recommended for the treatment of Inflammatory bowel disease, especially ulcerative colitis and Crohn's disease.

It binds to $\alpha 4\beta 7$ integrin, responsible of the immune cells (lymphocytes) attraction in the digestive tract and limits the intestinal inflammation which appears during these pathologies.

During the treatment, some patients can develop antibodies against Vedolizumab.

LISA-TRACKER Duo Vedolizumab () allows the detection of 2 parameters: Vedolizumab and anti-Vedolizumab antibodies. This kit allows the physician to monitor the level of these 2 parameters in patient sera.

SAMPLES COLLECTION AND HANDLING

- The test should be performed on plasma (heparin or citrate collection tubes) or on serum.
- Do not use plasma with EDTA collection tubes.
- Lipemic sera should be avoided, as well as samples which have been frozen and defrosted more than once.
- To avoid any non-specific binding, samples which have been frozen for more than 6 months or which are cloudy, should be centrifuged and filtered.

ASSAY PRINCIPLE

A. Dosage of Vedolizumab

The $\alpha 4\beta 7$ integrin is coated onto a polystyrene microtiter plate (6 strips of 8 wells).

- First, the diluted sample is added to the $\alpha 4\beta 7$ integrin coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Biotinylated antibody anti-vedolizumab is added. After incubation, unbound antibodies are removed by washing.
- Then horseradish peroxidase labelled streptavidin is added. The streptavidin binds to the complex formed with biotinylated anti-vedolizumab antibodies.

After incubation, the wells are washed again to eliminate any excess of conjugate.

- The bound enzyme is revealed by addition of substrate TMB (3,3',5,5' tetramethylbenzidin). The colour intensity is proportional to the amount of Vedolizumab.
- Adding H_2SO_4 allows to stop the enzymatic reaction.
- After stopping the reaction by H_2SO_4 , the optical density is read by a spectrophotometer at 450nm.

A range of calibration allows to define the quantity of Vedolizumab of each patient samples expressed in $\mu g/mL$.

B. Dosage of anti-Vedolizumab

The Vedolizumab is coated onto a polystyrene microtiter plate (6 strips of 8 wells).

- First, the diluted sample is added to the antibody coated well, which allows to bind. After incubation, unbound proteins are removed by washing.
- Biotinylated Vedolizumab is added. After incubation, unbound antibodies are removed by washing
- Then horseradish peroxidase labelled streptavidin is added. The streptavidin binds to the complex formed with biotinylated Vedolizumab. After incubation, the wells are washed again to eliminate any excess of conjugate.
- The bound enzyme is revealed by addition of substrate TMB (3,3',5,5' tetramethylbenzidin). The colour intensity is proportional to the amount of anti-Vedolizumab antibodies.
- Adding H_2SO_4 allows to stop the enzymatic reaction.
- After stopping the reaction by H_2SO_4 , the optical density is read by a spectrophotometer at 450nm.

A range of calibration allows to define the quantity of anti-Vedolizumab antibodies of each patient samples expressed in ng/mL .

REAGENTS

3 reagent families :

Color	Vedolizumab reagents	anti-Vedolizumab antibodies reagents	Common reagents
cap of vials	Blue	Yellow	Green, White, Black or Purple
microwells	Blue	Yellow	-

A) Specific reagents for the Vedolizumab determination

Strips of individual breakaway blue wells coated with human $\alpha 4\beta 7$ integrin. MP	6 strips
5 vials of « Vedolizumab » Standards, ($\mu\text{g/mL}$). <u>Ready to use.</u> The vials can be reused several times. The quantity of Vedolizumab is indicated on the vial label. Blue caps VED CAL n	5 x 1.5mL
«Positive control - Vedolizumab», ($\mu\text{g/mL}$). <u>To dilute.</u> VED CONTROL + The vials can be reused several times. The quantity of Vedolizumab is indicated on the vial label. Blue cap	1 x 250 μL
Biotinylated antibody vial. <u>Ready to use.</u> VED Ab BIOT Blue cap	1 x 7.5mL

B) Specific reagents for the anti-Vedolizumab antibodies determination

Strips of individual breakaway yellow wells coated with Vedolizumab. MP	6 strips
5 vials of « anti-Vedolizumab » Standards, (ng/mL). <u>Ready to use.</u> The vials can be reused several times. The quantity of anti-Vedolizumab is indicated on the vial label. Yellow caps A-VED CAL n	5 x 1.5mL
« Positive control – anti-Vedolizumab», (ng/mL). <u>To dilute.</u> A-VED CONTROL + The vial can be reused several times. The quantity of anti-Vedolizumab is indicated on the vial label. Yellow cap	1 x 1mL
Biotinylated antibody vial. <u>Ready to use.</u> A-VED Ab BIOT Yellow cap	1 x 7.5mL

C) Common reagents

HRP labelled Streptavidin. <u>Ready to use.</u> CONJ HRP Green cap	1 x 12mL
Phosphate-Tween Buffer pH 7,2 (10x) – <u>To reconstitute with distilled water.</u> White cap BUF WASH 10x	1 x 100mL
Substrate (TMB). <u>Ready to use.</u> SUBS TMB Black cap	1 x 12mL
Stop solution - H_2SO_4 (0.25 N). <u>Ready to use.</u> SOLN STOP Purple cap	1 x 15mL

MATERIAL REQUIRED BUT NOT PROVIDED

- distilled water
- precision pipettes
- microplate spectrophotometer with 450 nm filter
- 8 channel pipettes

STABILITY AND STORAGE

- Store reagents and micro-wells at $+2^\circ\text{C}/+8^\circ\text{C}$ in their own package.
- Do not use kits beyond the expiration date.
- Store unused strips into their plastic bag (closed securely) with the desiccant.
- Store all components immediately after use again at $+2^\circ\text{C}/+8^\circ\text{C}$.

SETUP

Except the TDL, which can be prepared in advance, all reagents must be prepared extemporaneously.

1. Dilution and Wash buffer (TDL)

- Dilute concentrated Phosphate-Tween Buffer 1/10 in distilled water.

BUF	WASH	10x
------------	-------------	------------

- Shelf life : 3 months at $+2^\circ\text{C}/+8^\circ\text{C}$ (avoid to use if signs of contamination or other visible changes occur).

NB. If there are crystals in the concentrated solution, warm the bottle up to $+37^\circ\text{C}$ for 15 minutes before use.

2. Preparation of samples and positive controls

a. Samples

- Vedolizumab determination

- Dilute to 1/401 in TDL
- Ex : 5 μL sample + 2mL TDL

- Vortex vigorously.

- Anti-Vedolizumab determination

- Dilute to 1/2 in TDL
- Ex : 130 μL sample + 130 μL TDL
- Vortex vigorously.

b. Positive controls

- Vedolizumab determination

- Dilute to 1/401 in TDL
- Ex : 5 μL positive control + 2mL TDL
- Vortex vigorously.

- anti-Vedolizumab determination

- Dilute to 1/2 in TDL
- Ex : 130 μL positive control + 130 μL TDL
- Vortex vigorously.

3. Use of ready-to-use biotinylated antibody.

- Estimate the amount required for handling and transfer to a tube.

4. Use of ready-to-use HRP Streptavidin conjugate.

- Estimate the amount required for handling and transfer to a tube.

5. Use of ready-to-use substrate.

- Estimate the amount required for handling and transfer to a dark tube.

PRECAUTIONS

Unpack all reagents in order to let them warm at room temperature (+18°C/+25°C) at least half an hour before starting the test.

⚠ The temperature of the reagents can impact the final result.

Check that all plates are well drained after each wash.

Avoid to use reagents if signs of contamination or other visible changes occur.

Human sources for the preparation of standards and controls have been tested and found negative for antibody to HIV 1 and 2, antibody to hepatitis C virus and hepatitis B virus antigen. Nevertheless, no test can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent. Therefore, the reagents should be handled as potentially infective materials.

Reagents in solution (except for substrate buffer and stop solution) contain <0.1% of sodium azide and <0.6% of ProClin® 300. Do not eat and avoid contact with skin and eyes. Azide can form explosive mixtures in copper or lead piping. Rinse thoroughly after flushing.

-  At this concentration, ProClin® 300 is irritating to eyes and skin, and may be detrimental if enough quantity is ingested. It is a skin sensitizer; prolonged or repeated exposure may cause allergic reaction in certain sensitive individuals.

LISA-TRACKER Duo Vedolizumab () has been developed according to EC regulation n°1272/2008 relating to the classification, labelling and packaging of substances and mixtures.

LISA-TRACKER Duo Vedolizumab () has been optimized for the use as described in this procedure. Do not substitute other manufacturer's reagents. Dilution or adulteration of these reagents may also affect the performance of the test. Close adherence to the test procedure will assure optimal performance.

METHOD

1. Preparing the test

Use the work sheet to note the sample locations.

Set out:

- 5 "standard" wells
- 1 well for positive control
- 1 well for each sample

For a simultaneous testing of the 2 parameters, unlink the exact number of wells needed. Return unused wells to plastic pouch provided in the kit, with the desiccant bag.

Remark :

If a dispensing/diluting device is used, place the specific wells in the following order : dosage of Vedolizumab then dosage of anti-Vedolizumab.

2. Samples, positive controls and standards incubation

Add 100 µL of standards, diluted controls or samples.

Incubate for 60 minutes at room temperature (+18°C/+25°C).

Wash step:

Remove the content of the wells by rapid inversion.

Wash 3 times with 300µL of TDL buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

3. Incubation of biotinylated antibodies

Add 100µL of **specific** biotinylated antibodies in identified wells.

Incubate for 60 minutes at room temperature (+18°C/+25°C).

Wash step:

Remove the content of the wells by rapid inversion.

Wash 3 times with 300µL of dilution and washing buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

4. Incubation of Conjugate

Add 100µL of conjugate.

Incubate for 30 minutes at room temperature (+18°C/+25°C).

Wash step:

Remove the content of the wells by rapid inversion.

Wash 3 times with 300µL of dilution and washing buffer.

Dry the microplate by tapping it gently on an absorbent paper to eliminate the excess of liquid.

5. Incubation of Substrate

Add 100µL substrate into each well.

Incubate for 15 minutes at room temperature (+18°C/+25°C), in the dark.

6. Stop of the reaction

Add 100µL of H₂SO₄ to each well.

7. Reading

Read the optical density of each well at 450nm within 30 minutes after stopping reaction.

RESULTS AND INTERPRETATION

A. Dosage of Vedolizumab

- The OD of the standard 1 should be at least 0.8.
- The positive control value should be comprised into the range indicated on the vial label.
- Trace a *degree 4 polynomial standard curve or a 4PL standard curve*, with plotting the units of the 5 standard points (µg/mL) along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).
- The Vedolizumab value can be directly read on the curve.
- Samples with values greater than that of standard 1 may be diluted to obtain a more precise result. The number of units should be multiplied by the selected dilution.

B. Dosage of anti-Vedolizumab

- The OD of the standard 1 should be at least 0.8.
- The positive control value should be comprised into the range indicated on the vial label.
- Trace the standard curve (polynomial curve), plotting the units of the 5 standard points (ng/mL) along the abscissa (X axis) and the corresponding OD values along the ordinate (Y axis).
- The anti-Vedolizumab value can be directly read on the curve.
- Samples with values greater than that of standard 1 may be diluted to obtain a more precise result. The number of units should be multiplied by the selected dilution.

CHARACTERISTICS AND PERFORMANCE OF THE TEST

Limits of detection / threshold values

Estimated on healthy patient population:

- 149 samples for Vedolizumab.
- 152 samples for anti-Vedolizumab.

Limit of detection Vedolizumab	Limit of detection anti-Vedolizumab
2 µg/mL 98 th percentile	35 ng/mL >99 th percentile

Assay range

Vedolizumab	anti-Vedolizumab
2 µg/mL – 60 µg/mL	35 ng/ml – 500 ng/ml

Interfering Substances Study

LISA-TRACKER Duo Vedolizumab () was evaluated to assess potential cross reactivity to other antibodies and interference from serum components (cryoglobulins, rheumatoid factors, heterophilic antibodies, high levels of triglycerides, bilirubin, IgG and/or IgM, and C1q proteins, autoantibodies) and on a population of patients' serums presenting inflammatory bowel disease. No interference was detected.

Precision

Parameters	Intra-run (30 tests in a same assay)		Inter-runs (20 different assays for parameter Vedolizumab) (11 different assays for parameter anti-Vedolizumab)	
	Mean	CV (%)	Mean	CV (%)
Vedolizumab (µg/ml)	8.4	8.1 %	4.3	15.2 %
	32.9	16.6 %	29.4	13.2 %
	54.5	7.5 %	46.5	12.4 %
anti-Vedolizumab (ng/ml)	54	15.7 %	61	14.7 %
	230	14.2 %	244	16.2 %
	399	10.5 %	360	11.7 %

LIMITS

The presence of biotin in patients' specimens can potentially impact the immunodosages using the Streptavidin-Biotin technology.

QUALITY CONTROL

It is recommended to use internally and externally sourced control material for the different specificities. **IMMUNO-TROL Vedolizumab** or **anti-Vedolizumab** controls () Cat.n: LTV 002-PC or LTV 003-PC, contain **Vedolizumab** or antibodies directed against **Vedolizumab**. These materials are to be assayed in the same manner as the unknown samples.

REFERENCES

Cherry LN & al., *Vedolizumab: an α4β7 integrin antagonist for ulcerative colitis and Crohn's disease*, Therapeutic Advances in Chronic Disease, Vol.6(5), 224-233, 2015.

Feagan BG & al., *Induction and Maintenance Therapy with Vedolizumab, a novel biologic therapy for ulcerative colitis*, Gastroenterology & Hepatology, 10, Issue 1, 2014.

Feagan BG & al., *Vedolizumab as induction and maintenance therapy for ulcerative colitis*, The new England journal of medicine, 369, 699-710, 2013.

McLean LP & al., *Vedolizumab for the treatment of ulcerative colitis and Crohn's disease*, Immunotherapy, 4, 883-898, 2012.

Petrovic A. & al., *LPAM (α4β7 integrin) is an important homing integrin on alloreactive T cells in the development of intestinal graft-versus-host disease*, Blood, 103, 4, 2004.

Piketty ML & al. Biotin: an emerging analytical interference, Ann Biol Clin2017 ; 75(4) : 366-8.

Rosario M. & al., *Population pharmacokinetics-pharmacodynamics of vedolizumab in patients with ulcerative colitis and Crohn's disease*, Alimentary Pharmacology and Therapeutics, 2015.

Sandborn WJ. & al., *Vedolizumab as induction and maintenance therapy for crohn's disease*, The new England journal of medicine, 369, 711-721, 2013.

Soler D. & al., *The binding specificity and selective antagonism of vedolizumab, an anti- α4β7integrin therapeutic antibody in development for inflammatory bowel diseases*,The journal of pharmacology and experimental therapeutics, 330, 864-875, 2009.

Yu Jin & al., *Meta-analysis of the effectiveness and safety of vedolizumab for ulcerative colitis*, World Journal of Gastroenterology, 21, issue 20, 2015.

SUMMARY OF METHOD

A) Sample Dilution

Vedolizumab	anti-Vedolizumab
1/401	1/2

B) Positive Control Dilution

Vedolizumab	anti-Vedolizumab
1/401	1/2

C) Procedure

Reagents	Procedure
Standards	100µL / wells
Diluted positive controls	
Diluted samples	
Incubation	1 h at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300µL / wells
Biotinylated antibodies	100µL / wells (specific reagents)
Incubation	1 h at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300µL / wells
HRP-Streptavidin	100µL / wells
Incubation	30 minutes at room temperature
Washing*	Wash 3 times with TDL buffer : 3 x 300µL / wells
Substrate (TMB)	100µL / wells
Incubation	15 minutes at room temperature. Protect from light.
Stop solution	100µL / wells

* Dry the microplate by tapping it gently on a towel to eliminate the excess of liquid.

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Vedolizumab assay

anti-Vedolizumab assay

c. 26 sera for anti-Vedolizumab determination

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